

AWP1 (m): 293T Lysate: sc-126472

BACKGROUND

AWP1 (associated with PKN protein), also known as zinc finger A20 domain-containing protein 3 (ZA20D3) or AN1-type zinc finger protein 6 (ZFAND6 or ZFAND5B), is a ubiquitously expressed protein with highest expression levels in liver, kidney, placenta, heart and skeletal muscle. AWP1 contains two conserved domains (one AN1-type zinc finger and one A20-type zinc finger), two PEST sequences, seven casein kinase II phosphorylation sites, two N-myristoylation sites and four N-glycosylation sites. The AN1-type zinc finger and four cysteine residues near the N-terminus are conserved between AWP1, OSISAP1 and ZNF216. Human AWP1 shares 55% homology with ZNF216 and mouse AWP1. It is a potent inhibitor of NFκB and it interacts with PKN, possibly participating in the regulation of PKN signal transduction pathways. Two AWP1 isoforms exist, due to alternative splicing. Isoform 1 is the mature full length protein and isoform 2 lacks amino acids 52-63.

REFERENCES

- Williams, K.J., Fisher, J.M. and Langridge, P. 1996. Development of a PCR-based allele-specific assay from an RFLP probe linked to resistance to cereal cyst nematode in wheat. *Genome* 39: 798-801.
- Duan, W., Sun, B., Li, T.W., Tan, B.J., Lee, M.K. and Teo, T.S. 2000. Cloning and characterization of AWP, a novel protein that associates with serine/threonine kinase PRK1 *in vivo*. *Gene* 256: 113-121.
- Mukhopadhyay, A., Vij, S. and Tyagi, A.K. 2004. Overexpression of a zinc-finger protein gene from rice confers tolerance to cold, dehydration, and salt stress in transgenic tobacco. *Proc. Natl. Acad. Sci. USA* 101: 6309-6314.
- Cao, Y.K., Mo, Y.Y., Tian, F.Z., Liu, Y.W., Deng, P., Qin, Q.H. and Jiang, Y. 2005. Construction of GFP-AWP1 fusion gene vector and its expression in 293 cells. *Di Yi Jun Yi Da Xue Xue Bao* 25: 174-176, 180.
- Lowes, D.A., Galley, H.F., Lowe, P.R., Rikke, B.A., Johnson, T.E. and Webster, N.R. 2005. A microarray analysis of potential genes underlying the neurosensitivity of mice to propofol. *Anesth. Analg.* 101: 697-704.
- Diatchenko, L., Romanov, S., Malinina, I., Clarke, J., Tchivilev, I., Li, X. and Makarov, S.S. 2005. Identification of novel mediators of NFκB through genome-wide survey of monocyte adherence-induced genes. *J. Leukoc. Biol.* 78: 1366-1377.
- Dash, D.P., Silvestri, G. and Hughes, A.E. 2006. Fine mapping of the keratoconus with cataract locus on chromosome 15q and candidate gene analysis. *Mol. Vis.* 12: 499-505.
- Wullaert, A., Heyninck, K., Janssens, S. and Beyaert, R. 2006. Ubiquitin: tool and target for intracellular NFκB inhibitors. *Trends Immunol.* 27: 533-540.
- Hishiya, A., Iemura, S., Natsume, T., Takayama, S., Ikeda, K. and Watanabe, K. 2006. A novel ubiquitin-binding protein ZNF216 functioning in muscle atrophy. *EMBO J.* 25: 554-564.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: Zfand6 (mouse) mapping to 7 D3.

PRODUCT

AWP1 (m): 293T Lysate represents a lysate of mouse AWP1 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

AWP1 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive AWP1 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.


AWP1 (B-5): sc-398513 is recommended as a positive control antibody for Western Blot analysis of enhanced mouse AWP1 expression in AWP1 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:

1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



AWP1 (B-5): sc-398513. Western blot analysis of AWP1 expression in non-transfected: sc-117752 (A) and mouse AWP1 transfected: sc-126472 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.