

Bag-1 (h2): 293T Lysate: sc-128079

BACKGROUND

The Bcl-2 family of proteins is characterized by its ability to modulate cell death (apoptosis) under a broad range of physiologic conditions. Bcl-2 and several related proteins function to inhibit apoptosis, while other members of the Bcl-2 family, such as Bax and Bak, enhance cell death under various conditions. For instance, Bcl-x_L represses cell death, while its shorter form, Bcl-x_S, promotes apoptosis. Dimerization of another member of this family, Bad, with Bcl-x_L, results in displacement of Bax from Bcl-x_L/Bax complexes and restoration of Bax-mediated apoptosis. A Bcl-2-binding protein, designated Bag-1, lacks significant homology with Bcl-2 or with other Bcl-2-related proteins. Bag-1 appears to function to enhance Bcl-2 protection from cell death, suggesting that Bag-1 represents a new type of anti-cell death gene. This also suggests that certain routes of apoptosis induction, previously ascribed to Bcl-2-independent pathways, may instead reflect a requirement for a combination of Bcl-2 and Bag-1.

REFERENCES

1. Nunez, G., London, L., Hockenbery, D., Alexander, M. and McKeown, J.P. 1990. Deregulated Bcl-2 gene expression selectively prolongs survival of growth factor-deprived hemopoietic cell lines. *J. Immunol.* 144: 3602-3610.
2. Hockenbery, D.M., Zutter, M., Hickey, W., Nahm, M. and Korsmeyer, S.J. 1991. Bcl-2 protein is topographically restricted in tissues characterized by apoptotic cell death. *Proc. Natl. Acad. Sci. USA* 88: 6961-6965.
3. Oltvai, Z.N., Millman, C.L. and Korsmeyer, S.J. 1993. Bcl-2 heterodimerizes *in vivo* with a conserved homolog, Bax, that accelerates programmed cell death. *Cell* 74: 609-619.
4. Yin, X.M., Oltvai, Z.N. and Korsmeyer, S.J. 1994. BH1 and BH2 domains of Bcl-2 are required for inhibition of apoptosis and heterodimerization with Bax. *Nature* 369: 321-323.
5. Chittenden, T., Harrington, E.A., O'Connor, R., Flemington, C., Lutz, R.J., Evan, G.I. and Guild, B.C. 1995. Induction of apoptosis by the Bcl-2 homologue Bak. *Nature* 374: 733-736.
6. Kiefer, M.C., Brauer, M.J., Powers, V.C., Wu, J.J., Umansky, S.R., Tomel, L.D. and Barr, P.J. 1995. Modulation of apoptosis by the widely distributed Bcl-2 homologue Bak. *Nature* 374: 736-739.
7. Takayama, S., Sato, T., Krajewski, S., Kochel, K., Irie, S., Millan, J.A. and Reed, J.C. 1995. Cloning and functional analysis of Bag-1: a novel Bcl-2-binding protein with anti-cell death activity. *Cell* 80: 279-284.
8. Yang, E., Zha, J., Jockel, J., Boise, L.H., Thompson, C.B. and Korsmeyer, S.J. 1995. Bad, a heterodimeric partner for Bcl-x_L and Bcl-2, displaces Bax and promotes cell death. *Cell* 80: 285-291.
9. Ishioka, S., Ezaka, Y., Umemura, K., Hayashi, T., Endo, T. and Saito, T. 2006. Proteomic analysis of mechanisms of hypoxia-induced apoptosis in trophoblastic cells. *Int. J. Med. Sci.* 4: 36-44.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: BAG1 (human) mapping to 9p13.3.

PRODUCT

Bag-1 (h2): 293T Lysate represents a lysate of human Bag-1 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

Bag-1 (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive Bag-1 antibodies. Recommended use: 10-20 µl per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

Bag-1 (SPM232): sc-56005 is recommended as a positive control antibody for Western Blot analysis of enhanced human Bag-1 expression in Bag-1 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).


RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



Bag-1 (SPM232): sc-56005. Western blot analysis of Bag-1 expression in non-transfected: sc-117752 (A) and human Bag-1 transfected: sc-128079 (B) 293T whole cell lysates.



Bag-1 (3.10G3E2): sc-56003. Western blot analysis of Bag-1 expression in non-transfected: sc-117752 (A) and human Bag-1 transfected: sc-128079 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.