

FBXO11 (LA-58): sc-130473

BACKGROUND

F-box proteins are critical components of the SCF (Skp1-CUL-1-F-box protein) type E3 ubiquitin ligase complex and are involved in substrate recognition and recruitment for ubiquitination. They are members of a larger family of proteins that are involved in the regulation of a wide variety of cellular processes (including the cell cycle, immune responses, signaling cascades and developmental events) through the targeting of proteins, such as cyclins, cyclin-dependent kinase inhibitors, I κ B- α and β -catenin, for proteasomal degradation. FBXO11 (F-box only protein 11), also known as VIT1 (vitiligo-associated protein 1), is a 927 amino acid nuclear protein that contains one UBR-type zinc finger, one F-box domain and 19 PbH1 repeats. Involved in protein ubiquitination, FBXO11 functions as a substrate recognition component of the SCF complex and is thought to bind to and inhibit the transcriptional activity of p53. Reduced expression of FBXO11 is associated with vitiligo, a disease characterized by progressive skin depigmentation. Multiple isoforms of FBXO11 exist due to alternative splicing events.

REFERENCES

- Cenciarelli, C., et al. 1999. Identification of a family of human F-box proteins. *Curr. Biol.* 9: 1177-1179.
- Winston, J.T., et al. 1999. A family of mammalian F-box proteins. *Curr. Biol.* 9: 1180-1182.
- Le Poole, I.C., et al. 2001. "VIT1", a novel gene associated with vitiligo. *Pigment Cell Res.* 14: 475-484.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607871. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Jin, J., et al. 2004. Systematic analysis and nomenclature of mammalian F-box proteins. *Genes Dev.* 18: 2573-2580.
- Segade, F., et al. 2006. Association of the FBXO11 gene with chronic otitis media with effusion and recurrent otitis media: the Minnesota COME/ROM family Study. *Arch. Otolaryngol. Head Neck Surg.* 132: 729-733.
- Cook, J.R., et al. 2006. FBXO11/PRMT9, a new protein arginine methyltransferase, symmetrically dimethylates arginine residues. *Biochem. Biophys. Res. Commun.* 342: 472-481.
- Abida, W.M., et al. 2007. FBXO11 promotes the Neddylation of p53 and inhibits its transcriptional activity. *J. Biol. Chem.* 282: 1797-1804.
- Lee, M.J., et al. 2008. Synthetic heterovalent inhibitors targeting recognition E3 components of the N-end rule pathway. *Proc. Natl. Acad. Sci. USA* 105: 100-105.

CHROMOSOMAL LOCATION

Genetic locus: FBXO11 (human) mapping to 2p16.3; Fbxo11 (mouse) mapping to 17 E4.

SOURCE

FBXO11 (LA-58) is a mouse monoclonal antibody raised against recombinant FBXO11 of human origin.

PRODUCT

Each vial contains 100 μ g IgG_{2a} kappa light chain in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

FBXO11 (LA-58) is recommended for detection of FBXO11 of mouse, rat and human origin by immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for FBXO11 siRNA (h): sc-94892, FBXO11 siRNA (m): sc-145102, FBXO11 shRNA Plasmid (h): sc-94892-SH, FBXO11 shRNA Plasmid (m): sc-145102-SH, FBXO11 shRNA (h) Lentiviral Particles: sc-94892-V and FBXO11 shRNA (m) Lentiviral Particles: sc-145102-V.


Molecular Weight of FBXO11: 103 kDa.

Molecular Weight of FBXO11 fragment: 14 kDa.

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850. 2) Immunohistochemistry: use m-IgG κ BP-HRP: sc-516102 with DAB, 50X: sc-24982 and Immunohistomount: sc-45086, or Organo/Limonene Mount: sc-45087.

DATA



FBXO11 (LA-58): sc-130473. Immunoperoxidase staining of formalin-fixed, paraffin-embedded human melanoma tissue showing nuclear and cytoplasmic localization.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.