

Nek5 (OC-65): sc-130492

BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including cell division, homeostasis and apoptosis. A group of proteins that are intimately involved in this process are the serine/threonine (Ser/Thr) protein kinases. Nek5 (NimA-related protein kinase 5) is a 708 amino acid protein that is related to NIMA, a protein that was originally discovered in *Aspergillus nidulans* and is necessary for entry into mitosis. One of several members of the Ser/Thr protein kinase super family, Nek5 contains one protein kinase domain through which it catalyzes the ATP-dependent phosphorylation of target proteins. Like NIMA, Nek5 may be involved in mitotic regulation and cell cycle control.

REFERENCES

- Hanks, S.K., Quinn, A.M. and Hunter, T. 1988. The protein kinase family: conserved features and deduced phylogeny of the catalytic domains. *Science* 241: 42-52.
- Lu, K.P. and Hunter, T. 1995. Evidence for a NIMA-like mitotic pathway in vertebrate cells. *Cell* 81: 413-424.
- Pu, R.T. and Osmani, S.A. 1995. Mitotic destruction of the cell cycle regulated NIMA protein kinase of *Aspergillus nidulans* is required for mitotic exit. *EMBO J.* 14: 995-1003.
- Lu, K.P. and Hunter, T. 1995. The NIMA kinase: a mitotic regulator in *Aspergillus nidulans* and vertebrate cells. *Prog. Cell Cycle Res.* 1: 187-205.
- Li, J.J. and Li, S.A. 2006. Mitotic kinases: the key to duplication, segregation, and cytokinesis errors, chromosomal instability, and oncogenesis. *Pharmacol. Ther.* 111: 974-984.
- O'regan, L., Blot, J. and Fry, A.M. 2007. Mitotic regulation by NIMA-related kinases. *Cell Div.* 2: 25.
- Vigneault, F., Lachance, D., Cloutier, M., Pelletier, G., Levasseur, C. and Séguin, A. 2007. Members of the plant NIMA-related kinases are involved in organ development and vascularization in poplar, *Arabidopsis* and rice. *Plant J.* 51: 575-588.
- Salaun, P., Rannou, Y. and Prigent, C. 2008. Cdk1, Plks, Auroras, and Neks: the mitotic bodyguards. *Adv. Exp. Med. Biol.* 617: 41-56.
- White, M.C. and Quarumby, L.M. 2008. The NIMA-family kinase, Nek1 affects the stability of centrosomes and ciliogenesis. *BMC Cell Biol.* 9: 29.

CHROMOSOMAL LOCATION

Genetic locus: NEK5 (human) mapping to 13q14.3; Nek5 (mouse) mapping to 8 A2.

SOURCE

Nek5 (OC-65) is a mouse monoclonal antibody raised against recombinant Nek5 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 100 µg IgG₁ kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

Nek5 (OC-65) is recommended for detection of Nek5 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)].

Suitable for use as control antibody for Nek5 siRNA (h): sc-75895, Nek5 siRNA (m): sc-149906, Nek5 shRNA Plasmid (h): sc-75895-SH, Nek5 shRNA Plasmid (m): sc-149906-SH, Nek5 shRNA (h) Lentiviral Particles: sc-75895-V and Nek5 shRNA (m) Lentiviral Particles: sc-149906-V.


Molecular Weight of Nek5: 81 kDa.

Positive Controls: IMR-32 cell lysate: sc-2409.

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



Nek5 (OC-65): sc-130492. Western blot analysis of Nek5 expression in IMR-32 whole cell lysate.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.