

DP-3 (8G11): sc-134321

BACKGROUND

DP-3, also known as TFDP3 (transcription factor Dp family, member 3), HCA661 (hepatocellular carcinoma-associated antigen 661) or CT30 (cancer/testis antigen 30), is a 405 amino acid protein that localizes to the nucleus and belongs to the E2F/DP family. Expressed predominately in testis, but also present at low levels in other tissues, DP-3 functions to regulate the activity of E2F transcription factor proteins and, via this regulation, plays a role in E2F-mediated apoptosis and may be involved in tumor growth and metastasis. The gene encoding DP-3 maps to human chromosome X, which contains nearly 153 million base pairs and houses over 1,000 genes. In conjunction with chromosome Y, chromosome X is responsible for sex determination. There are a number of conditions related to an abnormal number and combination of sex chromosomes, some of which include Turner's syndrome, color blindness, hemophilia and Duchenne muscular dystrophy.

REFERENCES

- Ginsberg, D. 2002. E2F1 pathways to apoptosis. *FEBS Lett.* 529: 122-125.
- Wang, Y., Han, K.J., Pang, X.W., Vaughan, H.A., Qu, W., Dong, X.Y., Peng, J.R., Zhao, H.T., Rui, J.A., Leng, X.S., Cebon, J., Burgess, A.W. and Chen, W.F. 2002. Large scale identification of human hepatocellular carcinoma-associated antigens by autoantibodies. *J. Immunol.* 169: 1102-1109.
- Bernardino-Sgherri, J., Flagiello, D. and Dutrillaux, B. 2002. Overall DNA methylation and chromatin structure of normal and abnormal X chromosomes. *Cytogenet. Genome Res.* 99: 85-91.
- Muntoni, F., Torelli, S. and Ferlini, A. 2003. Dystrophin and mutations: one gene, several proteins, multiple phenotypes. *Lancet Neurol.* 2: 731-740.
- Hitchens, M.R. and Robbins, P.D. 2003. The role of the transcription factor DP in apoptosis. *Apoptosis.* 8: 461-468.
- Deeb, S.S. 2005. The molecular basis of variation in human color vision. *Clin. Genet.* 67: 369-377.
- Hayashi, T., Kubo, A., Takeuchi, T., Gekka, T., Goto-Omoto, S. and Kitahara, K. 2006. Novel form of a single X-linked visual pigment gene in a unique dichromatic color-vision defect. *Vis. Neurosci.* 23: 411-417.
- Tian, C., Lv, D., Qiao, H., Zhang, J., Yin, Y.H., Qian, X.P., Wang, Y.P., Zhang, Y. and Chen, W.F. 2007. TFDP3 inhibits E2F1-induced, p53-mediated apoptosis. *Biochem. Biophys. Res. Commun.* 361: 20-25.
- Qiao, H., Di Stefano, L., Tian, C., Li, Y.Y., Yin, Y.H., Qian, X.P., Pang, X.W., Li, Y., McNutt, M.A., Helin, K., Zhang, Y. and Chen, W.F. 2007. Human TFDP3, a novel DP protein, inhibits DNA binding and transactivation by E2F. *J. Biol. Chem.* 282: 454-466.

CHROMOSOMAL LOCATION

Genetic locus: TFDP3 (human) mapping to Xq26.2.

SOURCE

DP-3 (8G11) is a mouse monoclonal antibody raised against recombinant DP-3 protein of human origin.

PRODUCT

Each vial contains 100 µg IgG_{2a} kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

DP-3 (8G11) is recommended for detection of DP-3 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for DP-3 siRNA (h): sc-106735, DP-3 shRNA Plasmid (h): sc-106735-SH and DP-3 shRNA (h) Lentiviral Particles: sc-106735-V.


Molecular Weight of DP-3: 45 kDa.

Positive Controls: IMR-32 cell lysate: sc-2409.

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



DP-3 (8G11): sc-134321. Western blot analysis of DP-3 expression in IMR-32 whole cell lysate.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.