

# L-FABP (2G4): sc-135628

## BACKGROUND

Fatty acid-binding proteins, designated FABPs, are a family of homologous cytoplasmic proteins that are expressed in a highly tissue-specific manner and play an integral role in the balance between lipid and carbohydrate metabolism. FABPs mediate fatty acid (FA) and/or hydrophobic ligand uptake, transport and targeting within their respective tissues. The mechanisms underlying these actions can give rise to both passive diffusional uptake and protein-mediated transmembrane transport of FAs. FABPs are expressed in adipocytes (A-FABP), brain (B-FABP), epidermis (E-FABP, also designated psoriasis-associated FABP or PA-FABP), muscle and heart (H-FABP, also designated mammary-derived growth inhibitor or MDGI), intestine (I-FABP), liver (L-FABP), myelin (M-FABP) and testis (T-FABP). Liver-specific FABP (L-FABP) expression is modulated by developmental, hormonal, dietary and pharmacological factors and is required for cholesterol synthesis and metabolism.

## REFERENCES

1. Kaikaus, R.M., Chan, W.K., Ortiz de Montellano, P.R. and Bass, N.M. 1993. Mechanisms of regulation of liver fatty acid-binding protein. *Mol. Cell. Biochem.* 123: 93-100.
2. Veerkamp, J.H. and Maatman, R.G. 1995. Cytoplasmic fatty acid-binding proteins: their structure and genes. *Prog. Lipid Res.* 34: 17-52.
3. Hotamisligil, G.S., Johnson, R.S., Distel, R.J., Ellis, R., Papaioannou, V.E. and Spiegelman, B.M. 1996. Uncoupling of obesity from Insulin resistance through a targeted mutation in aP2, the adipocyte fatty acid binding protein. *Science* 274: 1377-1379.
4. Murphy, E.J. 1998. L-FABP and I-FABP expression increase NBD-stearate uptake and cytoplasmic diffusion in L cells. *Am. J. Physiol.* 275: G244-G249.
5. Storch, J. and Thumser, A.E. 2000. The fatty acid transport function of fatty acid-binding proteins. *Biochim. Biophys. Acta* 1486: 28-44.
6. Glatz, J.F. and Storch, J. 2001. Unravelling the significance of cellular fatty acid-binding proteins. *Curr. Opin. Lipidol.* 12: 267-274.
7. Veerkamp, J.H. and Zimmerman, A.W. 2001. Fatty acid-binding proteins of nervous tissue. *J. Mol. Neurosci.* 16: 133-142.
8. Online Mendelian Inheritance in Man, OMIM™. 2001. Johns Hopkins University, Baltimore, MD. MIM Number: 134640. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

## CHROMOSOMAL LOCATION

Genetic locus: FABP1 (human) mapping to 2p11.2.

## SOURCE

L-FABP (2G4) is a mouse monoclonal antibody raised against full-length recombinant L-FABP of human origin.

## PRODUCT

Each vial contains 100 µg IgG<sub>1</sub> in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS


L-FABP (2G4) is recommended for detection of L-FABP of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for L-FABP siRNA (h): sc-41243, L-FABP shRNA Plasmid (h): sc-41243-SH and L-FABP shRNA (h) Lentiviral Particles: sc-41243-V.

Molecular Weight of L-FABP: 14 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227.

## DATA



L-FABP (2G4): sc-135628. Western blot analysis of L-FABP expression in Hep G2 whole cell lysate.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.