

GSKIP siRNA (m): sc-140326

BACKGROUND

GSKIP (GSK3-beta interaction protein), also known as C14orf129 or HSPC210, is a 139 amino acid protein belonging to the UPF0279 family. Localizing to the cytoplasm, GSKIP is expressed in heart, brain, placenta, liver, skeletal muscle, kidney, testis, lung and pancreas. GSKIP interacts directly with GSK-3B, a protein that plays an important role in various physiological functions and regulates axons and dendrites, resulting in GSK-3B inhibition. The gene encoding GSKIP maps to human chromosome 14, which contains approximately 700 genes and 106 million base pairs, and makes up about 3.5% of human cellular DNA. Chromosome 14 encodes the presenilin 1 (PSEN1) gene, which is one of the three key genes associated with the development of Alzheimer's disease. The genetic disorder α 1-antitrypsin deficiency is also associated with chromosome 14.

REFERENCES

1. Heilig, R., Eckenberg, R., Petit, J.L., Fonknechten, N., Da Silva, C., Cattolico, L., Levy, M., Barbe, V., de Berardinis, V., Ureta-Vidal, A., Pelletier, E., Vico, V., Anthouard, V., Rowen, L., Madan, A., Qin, S., Sun, H., et al. 2003. The DNA sequence and analysis of human chromosome 14. *Nature* 421: 601-607.
2. Avramopoulos, D., Fallin, M.D. and Bassett, S.S. 2005. Linkage to chromosome 14q in Alzheimer's disease (AD) patients without psychotic symptoms. *Am. J. Med. Genet. B Neuropsychiatr. Genet.* 132B: 9-13.
3. Chou, H.Y., Howng, S.L., Cheng, T.S., Hsiao, Y.L., Lieu, A.S., Loh, J.K., Hwang, S.L., Lin, C.C., Hsu, C.M., Wang, C., Lee, C.I., Lu, P.J., Chou, C.K., Huang, C.Y. and Hong, Y.R. 2006. GSKIP is homologous to the Axin GSK3 β interaction domain and functions as a negative regulator of GSK3 β . *Biochemistry* 45: 11379-11389.
4. Lerner, A.J. and Doran, M. 2009. Genotype-phenotype relationships of presenilin-1 mutations in Alzheimer's disease: an update. *J. Alzheimers Dis.* 17: 259-265.
5. Lin, C.C., Chou, C.H., Howng, S.L., Hsu, C.Y., Hwang, C.C., Wang, C., Hsu, C.M. and Hong, Y.R. 2009. GSKIP, an inhibitor of GSK3 β , mediates the N-cadherin/ β -catenin pool in the differentiation of SH-SY5Y cells. *J. Cell. Biochem.* 108: 1325-1336.
6. Topic, A., Alempijevic, T., Milutinovic, A.S. and Kovacevic, N. 2009. α 1-antitrypsin phenotypes in adult liver disease patients. *Ups. J. Med. Sci.* 114: 228-234.
7. Howng, S.L., Hwang, C.C., Hsu, C.Y., Hsu, M.Y., Teng, C.Y., Chou, C.H., Lee, M.F., Wu, C.H., Chiou, S.J., Lieu, A.S., Loh, J.K., Yang, C.N., Lin, C.S. and Hong, Y.R. 2010. Involvement of the residues of GSKIP, AxinGID, and FRATtide in their binding with GSK3beta to unravel a novel C-terminal scaffold-binding region. *Mol. Cell. Biochem.* 339: 23-33.
8. Tang, X.N., Lo, C.W., Chuang, Y.C., Chen, C.T., Sun, Y.C., Hong, Y.R. and Yang, C.N. 2011. Prediction of the binding mode between GSK3 β and a peptide derived from GSKIP using molecular dynamics simulation. *Biopolymers* 95: 461-471.

CHROMOSOMAL LOCATION

Genetic locus: Gskip (mouse) mapping to 12 E.

PRODUCT

GSKIP siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GSKIP shRNA Plasmid (m): sc-140326-SH and GSKIP shRNA (m) Lentiviral Particles: sc-140326-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GSKIP siRNA (m) is recommended for the inhibition of GSKIP expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GSKIP gene expression knockdown using RT-PCR Primer: GSKIP (m)-PR: sc-140326-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.