

α -protein kinase 1 siRNA (m): sc-140597

BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including cell division, homeostasis and apoptosis. α -protein kinase 1, also known as ALPK1, LAK or KIAA1527, is a 1,244 amino acid protein that contains one α -type protein kinase domain and belongs to the protein kinase superfamily. Expressed at high levels in liver, α -protein kinase 1 recognizes and phosphorylates specific phosphorylation sites that are surrounded by peptides which have an α -helical conformation. One of the targets of α -protein kinase 1 is Myosin I and, via its ability to phosphorylate Myosin I, α -protein kinase 1 is thought to play an important role in the apical trafficking of vesicles carrying raft-associated sucrose-isomaltase (SI).

REFERENCES

- Ryazanov, A.G., Pavur, K.S. and Dorovkov, M.V. 1999. α -kinases: a new class of protein kinases with a novel catalytic domain. *Curr. Biol.* 9: R43-R45.
- Nakajima, D., Okazaki, N., Yamakawa, H., Kikuno, R., Ohara, O. and Nagase, T. 2002. Construction of expression-ready cDNA clones for KIAA genes: manual curation of 330 KIAA cDNA clones. *DNA Res.* 9: 99-106.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607347. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Yamada, S., Ohira, M., Horie, H., Ando, K., Takayasu, H., Suzuki, Y., Sugano, S., Hirata, T., Goto, T., Matsunaga, T., Hiyama, E., Hayashi, Y., Ando, H., Suita, S., Kaneko, M., Sasaki, F., Hashizume, K., et al. 2004. Expression profiling and differential screening between hepatoblastomas and the corresponding normal livers: identification of high expression of the PLK1 oncogene as a poor-prognostic indicator of hepatoblastomas. *Oncogene* 23: 5901-5911.
- Drennan, D. and Ryazanov, A.G. 2004. α -kinases: analysis of the family and comparison with conventional protein kinases. *Prog. Biophys. Mol. Biol.* 85: 1-32.
- Heine, M., Cramm-Behrens, C.I., Ansari, A., Chu, H.P., Ryazanov, A.G., Naim, H.Y. and Jacob, R. 2005. α -kinase 1, a new component in apical protein transport. *J. Biol. Chem.* 280: 25637-25643.

CHROMOSOMAL LOCATION

Genetic locus: Alpk1 (mouse) mapping to 3 G2.

PRODUCT

protein kinase 1 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see α -protein kinase 1 shRNA Plasmid (m): sc-140597-SH and α -protein kinase 1 shRNA (m) Lentiviral Particles: sc-140597-V as alternate gene silencing products.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

α -protein kinase 1 siRNA (m) is recommended for the inhibition of α -protein kinase 1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor α -protein kinase 1 gene expression knockdown using RT-PCR Primer: α -protein kinase 1 (m)-PR: sc-140597-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.