



# ACN9 siRNA (m): sc-140810

## BACKGROUND

ACN9 (CN9 homolog (*S. cerevisiae*)), also known as protein ACN9 homolog, mitochondrial, is a 125 amino acid protein belonging to the ACN9 family. ACN9 homologs exist in human, mouse, and nematode. Encoded by a gene that maps to human chromosome 7q21.3, ACN9 participates in gluconeogenesis regulation and ethanol or acetate assimilation into carbohydrate. ACN9 also functions in activities related to the overexpression of oxidative enzymes, such as glyoxylate cycle activity and acetyl-CoA metabolism. Localizing to the mitochondrial intermembrane space, ACN9 expression can be slightly repressed by glucose. ACN9 may play a role in predisposition to alcohol dependence and may be linked to early-onset breast cancer.

## REFERENCES

- Whitford, M., Stewart, S., Kuzio, J. and Faulkner, P. 1989. Identification and sequence analysis of a gene encoding gp67, an abundant envelope glycoprotein of the baculovirus *Autographa californica* nuclear polyhedrosis virus. *J. Virol.* 63: 1393-1399.
- McCammon, M.T. 1996. Mutants of *Saccharomyces cerevisiae* with defects in acetate metabolism: isolation and characterization of Acn-mutants. *Genetics* 144: 57-69.
- Collier, L.S., Gaines, G.L. and Neidle, E.L. 1998. Regulation of benzoate degradation in *Acinetobacter* sp. strain ADP1 by BenM, a LysR-type transcriptional activator. *J. Bacteriol.* 180: 2493-2501.
- Dennis, R.A. and McCammon, M.T. 1999. ACN9 is a novel protein of gluconeogenesis that is located in the mitochondrial intermembrane space. *Eur. J. Biochem.* 261: 236-243.
- Nishizawa, M., Katou, Y., Shirahige, K. and Toh-e, A. 2004. Yeast Pho85 kinase is required for proper gene expression during the diauxic shift. *Yeast* 21: 903-918.
- Reinders, J., Zahedi, R.P., Pfanner, N., Meisinger, C. and Sickmann, A. 2006. Toward the complete yeast mitochondrial proteome: multidimensional separation techniques for mitochondrial proteomics. *J. Proteome Res.* 5: 1543-1554.
- Dick, D.M., Aliev, F., Wang, J.C., Saccone, S., Hinrichs, A., Bertelsen, S., Budde, J., Saccone, N., Foroud, T., Nurnberger, J., Xuei, X., Conneally, P.M., Schuckit, M., Almasy, L., Crowe, R., Kuperman, S., Kramer, J., et al. 2008. A systematic single nucleotide polymorphism screen to fine-map alcohol dependence genes on chromosome 7 identifies association with a novel susceptibility gene ACN9. *Biol. Psychiatry* 63: 1047-1053.
- Kibriya, M.G., Jasmine, F., Argos, M., Andrulis, I.L., John, E.M., Chang-Claude, J. and Ahsan, H. 2009. A pilot genome-wide association study of early-onset breast cancer. *Breast Cancer Res. Treat.* 114: 463-477.

## CHROMOSOMAL LOCATION

Genetic locus: Acn9 (mouse) mapping to 6 A1.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

ACN9 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ACN9 shRNA Plasmid (m): sc-140810-SH and ACN9 shRNA (m) Lentiviral Particles: sc-140810-V as alternate gene silencing products.

For independent verification of ACN9 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-140810A and sc-140810B.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

ACN9 siRNA (m) is recommended for the inhibition of ACN9 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ACN9 gene expression knockdown using RT-PCR Primer: ACN9 (m)-PR: sc-140810-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.