

ADAMTS-L5 siRNA (m): sc-140871

BACKGROUND

ADAMTS (a disintegrin and metalloproteinase domain with thrombospondin type-1 modules) is a family of zinc-dependent proteases that are implicated in a variety of normal and pathological conditions, including arthritis and cancer. ADAMTS protein family members contain an amino-terminal propeptide domain, a metalloproteinase domain, a disintegrin-like domain and a carboxy-terminus that contains a varying number of thrombospondin type 1 (TSP-1) motifs. ADAMTS-L5 (ADAMTS-like protein 5), also known as thrombospondin type-1 domain-containing protein 6 or THSD6, is a 471 amino acid secreted protein. Encoded by a gene that maps to human chromosome 19p13.3, ADAMTS-L5 contains a single TSP-1 domain and participates in protein and zinc ion binding. Although strongly similar to members of the ADAMTS family, ADAMTS-L5 lacks the metalloprotease and disintegrin-like domains typical of the family.

REFERENCES

1. Saban, M.R., Hellmich, H.L., Turner, M., Nguyen, N.B., Vadigepalli, R., Dyer, D.W., Hurst, R.E., Centola, M. and Saban, R. 2006. The inflammatory and normal transcriptome of mouse bladder detrusor and mucosa. *BMC Physiol.* 6: 1.
2. Bocker, C.N., Vasilou, V. and Nebert, D.W. 2009. Evolutionary divergence and functions of the ADAM and ADAMTS gene families. *Hum. Genomics* 4: 43-55.
3. Wang, L.W., Leonhard-Melief, C., Haltiwanger, R.S. and Apte, S.S. 2009. Post-translational modification of thrombospondin type-1 repeats in ADAMTS-like 1/punctin-1 by C-mannosylation of tryptophan. *J. Biol. Chem.* 284: 30004-30015.
4. Apte, S.S. 2009. A disintegrin-like and metalloprotease (reprolysin-type) with thrombospondin type 1 motif (ADAMTS) superfamily: functions and mechanisms. *J. Biol. Chem.* 284: 31493-31497.
5. Okaty, B.W., Miller, M.N., Sugino, K., Hempel, C.M. and Nelson, S.B. 2009. Transcriptional and electrophysiological maturation of neocortical fast-spiking GABAergic interneurons. *J. Neurosci.* 29: 7040-7052.
6. Paruthiyil, S., Cvorovic, A., Zhao, X., Wu, Z., Sui, Y., Staub, R.E., Baggett, S., Herber, C.B., Griffin, C., Tagliaferri, M., Harris, H.A., Cohen, I., Bjeldanes, L.F., Speed, T.P., Schaufele, F. and Leitman, D.C. 2009. Drug and cell type-specific regulation of genes with different classes of estrogen receptor β -selective agonists. *PLoS ONE* 4: e6271.
7. Sollazzo, V., Palmieri, A., Pezzetti, F., Massari, L. and Carinci, F. 2010. Effects of pulsed electromagnetic fields on human osteoblastlike cells (MG-63): a pilot study. *Clin. Orthop. Relat. Res.* 468: 2260-2277.
8. SWISS-PROT/TrEMBL (Q6ZMM2). World Wide Web URL: <http://www.uniprot.org/uniprot/Q6ZMM2>

CHROMOSOMAL LOCATION

Genetic locus: Adamtsl5 (mouse) mapping to 10 C1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

ADAMTS-L5 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ADAMTS-L5 shRNA Plasmid (m): sc-140871-SH and ADAMTS-L5 shRNA (m) Lentiviral Particles: sc-140871-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ADAMTS-L5 siRNA (m) is recommended for the inhibition of ADAMTS-L5 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ADAMTS-L5 gene expression knockdown using RT-PCR Primer: ADAMTS-L5 (m)-PR: sc-140871-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.