



AEBP2 siRNA (m): sc-140891

BACKGROUND

Adipocyte enhancer-binding protein 2 (AEBP2), also known as zinc finger protein AEBP2, is an 817 amino acid member of the AEBP2 C₂H₂-type zinc-finger protein family. Localized to the nucleus, AEBP2 acts as a DNA-binding transcriptional repressor of the adipose P2 (aP2) gene. The aP2 gene, which encodes the adipose fatty acid-binding protein, plays a critical role in triglyceride metabolism during adipocyte differentiation. The AE-1 region in the proximal promoter region of the aP2 gene functions as either a positive or negative regulatory element. C/EBP α binds to the AE-1 sequence and functions as a transcriptional activator of aP2, whereas other proteins, such as AEBP2, bind to the region and repress gene expression. AEBP2 contains three C₂H₂-type zinc fingers and it has been shown that not all of the zinc fingers are involved in DNA binding. Three isoforms of AEBP2 exist as a result of alternative splicing events.

REFERENCES

1. He, G.P., Muise, A., Li, A.W. and Ro, H.S. 1995. A eukaryotic transcriptional repressor with carboxypeptidase activity. *Nature* 378: 92-96.
2. Sandell, L.J. 1996. Genes and gene regulation of extracellular matrix proteins: an introduction. *Connect. Tissue Res.* 35: 1-6.
3. He, G.P., Kim, S. and Ro, H.S. 1999. Cloning and characterization of a novel zinc finger transcriptional repressor. A direct role of the zinc finger motif in repression. *J. Biol. Chem.* 274: 14678-14684.
4. Cowherd, R.M., Lyle, R.E. and McGehee, R.E. 1999. Molecular regulation of adipocyte differentiation. *Semin. Cell Dev. Biol.* 10: 3-10.
5. Ro, H.S., Kim, S.W., Wu, D., Webber, C. and Nicholson, T.E. 2001. Gene structure and expression of the mouse adipocyte enhancer-binding protein. *Gene* 280: 123-133.
6. Cao, R., Wang, L., Wang, H., Xia, L., Erdjument-Bromage, H., Tempst, P., Jones, R.S. and Zhang, Y. 2002. Role of histone H3 lysine 27 methylation in Polycomb-group silencing. *Science* 298: 1039-1043.
7. Cao, R. and Zhang, Y. 2004. SUZ12 is required for both the histone methyltransferase activity and the silencing function of the EED-EZH2 complex. *Mol. Cell* 15: 57-67.
8. Lyons, P.J., Muise, A.M. and Ro, H.S. 2005. MAPK modulates the DNA binding of adipocyte enhancer-binding protein 1. *Biochemistry* 44: 926-931.
9. Lyons, P.J., Mattatall, N.R. and Ro, H.S. 2006. Modeling and functional analysis of AEBP1, a transcriptional repressor. *Proteins* 63: 1069-1083.

CHROMOSOMAL LOCATION

Genetic locus: Aebp2 (mouse) mapping to 6 G2.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

AEBP2 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see AEBP2 shRNA Plasmid (m): sc-140891-SH and AEBP2 shRNA (m) Lentiviral Particles: sc-140891-V as alternate gene silencing products.

For independent verification of AEBP2 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-140891A and sc-140891B.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

AEBP2 siRNA (m) is recommended for the inhibition of AEBP2 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor AEBP2 gene expression knockdown using RT-PCR Primer: AEBP2 (m)-PR: sc-140891-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.