

▶ AIM1L siRNA (m): sc-140967

BACKGROUND

AIM1L (absent in melanoma 1-like), also known as CRYBG2 (β/γ crystallin domain-containing protein 2), is a 616 amino acid protein that participates in sugar binding. AIM1L contains ten β/γ crystallin "Greek key" domains and one ricin B-type lectin domain. Associated with tumorigenesis, AIM1L is commonly down-regulated in colorectal amilial adenomatous polyposis (FAP) polyps, compared to normal colon tissue. AIM1L is encoded by a gene that maps to human chromosome 1p36.11. As the largest human chromosome, chromosome 1 makes up approximately 8% of the human genome and contains 260 million base pairs encoding 3,000 genes. Numerous diseases are linked to chromosome 1, notably the rare aging disease Hutchinson-Gilford progeria, which is associated with Lamin A. When defective, Lamin A can accumulate in nucleus, causing characteristic nuclear blebs. The MUTYH gene is located on chromosome 1 and is partially responsible for familial adenomatous polyposis. Stickler syndrome, Parkinson's disease, Gaucher disease and Usher syndrome are also associated with chromosome 1. Aberrations in chromosome 1 exist in a variety of cancers, including head and neck cancer, malignant melanoma and multiple myeloma.

REFERENCES

1. Teichmann, U., Ray, M.E., Ellison, J., Graham, C., Wistow, G., Meltzer, P.S., Trent, J.M. and Pavan, W.J. 1998. Cloning and tissue expression of the mouse ortholog of AIM1, a β/γ -crystallin superfamily member. *Mamm. Genome* 9: 715-720.
2. Weise, A., Starke, H., Mrasek, K., Claussen, U. and Liehr, T. 2005. New insights into the evolution of chromosome 1. *Cytogenet. Genome Res.* 108: 217-222.
3. Marzin, Y., Jamet, D., Douet-Guilbert, N., Morel, F., Le Bris, M.J., Morice, P., Abgrall, J.F., Berthou, C. and De Braekeleer, M. 2006. Chromosome 1 abnormalities in multiple myeloma. *Anticancer Res.* 26: 953-959.
4. Gregory, S.G., Barlow, K.F., McLay, K.E., Kaul, R., Swarbreck, D., Dunham, A., Scott, C.E., Howe, K.L., Woodfine, K., Spencer, C.C., Jones, M.C., Gillson, C., Searle, S., Zhou, Y., et al. 2006. The DNA sequence and biological annotation of human chromosome 1. *Nature* 441: 315-321.
5. McClintock, D., Gordon, L.B. and Djabali, K. 2006. Hutchinson-Gilford progeria mutant lamin A primarily targets human vascular cells as detected by an anti-Lamin A G608G antibody. *Proc. Natl. Acad. Sci. USA* 103: 2154-2159.
6. Bowden, N.A., Croft, A. and Scott, R.J. 2007. Gene expression profiling in familial adenomatous polyposis adenomas and desmoid disease. *Hered. Cancer Clin. Pract.* 5: 79-96.
7. Jelinsky, S.A., Choe, S.E., Crabtree, J.S., Cotreau, M.M., Wilson, E., Saraf, K., Dorner, A.J., Brown, E.L., Peano, B.J., Zhang, X., Winneker, R.C. and Harris, H.A. 2008. Molecular analysis of the vaginal response to estrogens in the ovariectomized rat and postmenopausal woman. *BMC Med. Genomics* 1: 27.
8. Kirpich, I.A., Gobejishvili, L.N., Homme, M.B., Waigel, S., Cave, M., Arteel, G., Barve, S.S., McClain, C.J. and Deaciuc, I.V. 2010. Integrated hepatic transcriptome and proteome analysis of mice with high-fat diet-induced nonalcoholic fatty liver disease. *J. Nutr. Biochem.* 22: 38-45.

CHROMOSOMAL LOCATION

Genetic locus: Aim1l (mouse) mapping to 4 D3.

PRODUCT

AIM1L siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see AIM1L shRNA Plasmid (m): sc-140967-SH and AIM1L shRNA (m) Lentiviral Particles: sc-140967-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

AIM1L siRNA (m) is recommended for the inhibition of AIM1L expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor AIM1L gene expression knockdown using RT-PCR Primer: AIM1L (m)-PR: sc-140967-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.