

# ANKMY2 siRNA (m): sc-141072

## BACKGROUND

ANKMY2 (ankyrin repeat and MYND domain containing 2) is a 441 amino acid protein that contains three ANK repeats and one MYND-type zinc finger. Encoded by a gene that maps to human chromosome 7p21.1, ANKMY2 is conserved in chimpanzee, canine, bovine, mouse, chicken, zebrafish, fruit fly, mosquito and *Caenorhabditis elegans*. Downregulation of ANKMY2, associated with frequent deletions of human chromosome 7p22.1, indicate that ANKMY2 may play a role in the pathogenesis of natural killer (NK)-cell malignancies. ANKMY2 is also upregulated by enforced expression of Hox11, which functions broadly to hinder hemopoiesis, diverts differentiation to an alternative fate and promotes pre-leukaemic states. ANKMY2 may also participate in zinc ion binding, cell adhesion, cellular morphogenesis, neurite outgrowth and human cortical functions.

## REFERENCES

1. Oldham, M.C., Horvath, S. and Geschwind, D.H. 2006. Conservation and evolution of gene coexpression networks in human and chimpanzee brains. *Proc. Natl. Acad. Sci. USA* 103: 17973-17978.
2. Dixon, D.N., Izon, D.J., Dagger, S., Callow, M.J., Taplin, R.H., Kees, U.R. and Greene, W.K. 2007. TLX1/HOX11 transcription factor inhibits differentiation and promotes a non-haemopoietic phenotype in murine bone marrow cells. *Br. J. Haematol.* 138: 54-67.
3. Hesson, L.B., Cooper, W.N. and Latif, F. 2007. Evaluation of the 3p21.3 tumour-suppressor gene cluster. *Oncogene* 26: 7283-7301.
4. Mahler, K.L., Fleming, J.L., Dworkin, A.M., Gladman, N., Cho, H.Y., Mao, J.H., Balmain, A. and Toland, A.E. 2008. Sequence divergence of *Mus spretus* and *Mus musculus* across a skin cancer susceptibility locus. *BMC Genomics* 9: 626.
5. Delplanque, J., Devos, D., Vuillaume, I., De Beccdelievre, A., Vangelder, E., Maurage, C.A., Dujardin, K., Destee, A. and Sablonniere, B. 2008. Slowly progressive spinocerebellar ataxia with extrapyramidal signs and mild cognitive impairment (SCA21). *Cerebellum* 7: 179-183.
6. Haider, N.B., Mollema, N., Gaule, M., Yuan, Y., Sachs, A.J., Nystuen, A.M., Naggert, J.K. and Nishina, P.M. 2009. Nr2e3-directed transcriptional regulation of genes involved in photoreceptor development and cell-type specific phototransduction. *Exp. Eye Res.* 89: 365-372.
7. Iqbal, J., Kucuk, C., Deleew, R.J., Srivastava, G., Tam, W., Geng, H., Klinkebiel, D., Christman, J.K., Patel, K., Cao, K., Shen, L., Dybkaer, K., Tsui, I.F., Ali, H., Shimizu, N., Au, W.Y., Lam, W.L. and Chan, W.C. 2009. Genomic analyses reveal global functional alterations that promote tumor growth and novel tumor suppressor genes in natural killer-cell malignancies. *Leukemia* 23: 1139-1151.
8. Ray, A., Liu, J., Ayoubi, P. and Pope, C. 2010. Dose-related gene expression changes in forebrain following acute, low-level chlorpyrifos exposure in neonatal rats. *Toxicol. Appl. Pharmacol.* 248: 144-155.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## CHROMOSOMAL LOCATION

Genetic locus: Ankmy2 (mouse) mapping to 12 A3.

## PRODUCT

ANKMY2 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ANKMY2 shRNA Plasmid (m): sc-141072-SH and ANKMY2 shRNA (m) Lentiviral Particles: sc-141072-V as alternate gene silencing products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

ANKMY2 siRNA (m) is recommended for the inhibition of ANKMY2 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ANKMY2 gene expression knockdown using RT-PCR Primer: ANKMY2 (m)-PR: sc-141072-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.