

ANKS4B siRNA (m): sc-141120

BACKGROUND

ANKS4B (ankyrin repeat and sterile alpha motif domain containing 4B), also known as Harp (harmonin-interacting ankyrin repeat-containing protein), is a 460 amino acid coiled-coil protein that interacts with the Harmonin PDZ1 domain in epithelial tissues by way of its C-terminus. Encoded by a gene that maps to human chromosome 16p12.2, ANKS4B contains two exons, three ANK repeats and one SAM (sterile α motif) domain. ANKS4B is expressed in kidney, liver, small intestine, colon, pancreas and purified pancreatic cells, as well as developing lung, kidney, salivary glands and cochlea. Conserved in chimpanzee, canine, bovine, mouse, chicken and zebrafish, ANKS4B shares significant homology with USH1G, which plays a major role in Usher syndrome. ANKS4B may also play a role in Usher syndrome and may facilitate signal transduction in epithelia.

REFERENCES

1. Weil, D., El-Amraoui, A., Masmoudi, S., Mustapha, M., Kikkawa, Y., Laine, S., Delmaghani, S., Adato, A., Nadifi, S., Zina, Z.B., Hamel, C., Gal, A., Ayadi, H., Yonekawa, H. and Petit, C. 2003. Usher syndrome type I G (USH1G) is caused by mutations in the gene encoding SANS, a protein that associates with the USH1C protein, harmonin. *Hum. Mol. Genet.* 12: 463-471.
2. Johnston, A.M., Naselli, G., Niwa, H., Brodnicki, T., Harrison, L.C. and Góñez, L.J. 2004. Harp (harmonin-interacting, ankyrin repeat-containing protein), a novel protein that interacts with harmonin in epithelial tissues. *Genes Cells* 9: 967-982.
3. Reiners, J. and Wolfrum, U. 2006. Molecular analysis of the supramolecular usher protein complex in the retina. Harmonin as the key protein of the Usher syndrome. *Adv. Exp. Med. Biol.* 572: 349-353.
4. Yokota, T., Mishra, M., Akatsu, H., Tani, Y., Miyauchi, T., Yamamoto, T., Kosaka, K., Nagai, Y., Sawada, T. and Heese, K. 2006. Brain site-specific gene expression analysis in Alzheimer's disease patients. *Eur. J. Clin. Invest.* 36: 820-830.
5. Kremer, H., van Wijk, E., Märker, T., Wolfrum, U. and Roepman, R. 2006. Usher syndrome: molecular links of pathogenesis, proteins and pathways. *Hum. Mol. Genet.* 15: R262-R270.
6. Dawson, P.A., Gardiner, B., Grimmond, S. and Markovich, D. 2006. Transcriptional profile reveals altered hepatic lipid and cholesterol metabolism in hyposulfatemic NaS1 null mice. *Physiol. Genomics* 26: 116-124.
7. Martin, B., Brenneman, R., Becker, K.G., Gucek, M., Cole, R.N. and Maudsley, S. 2008. iTRAQ analysis of complex proteome alterations in 3xTgAD Alzheimer's mice: understanding the interface between physiology and disease. *PLoS ONE* 3: e2750.

CHROMOSOMAL LOCATION

Genetic locus: Anks4b (mouse) mapping to 7 F2.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

ANKS4B siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ANKS4B shRNA Plasmid (m): sc-141120-SH and ANKS4B shRNA (m) Lentiviral Particles: sc-141120-V as alternate gene silencing products.

For independent verification of ANKS4B (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-141120A, sc-141120B and sc-141120C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ANKS4B siRNA (m) is recommended for the inhibition of ANKS4B expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ANKS4B gene expression knockdown using RT-PCR Primer: ANKS4B (m)-PR: sc-141120-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.