

apoL-VI siRNA (m): sc-141162

BACKGROUND

The apolipoprotein L gene family maps to a region on chromosome 22 and encodes six highly homologous proteins, designated apoL-I, apoL-II, apoL-III, apoL-IV, apoL-V and apoL-VI, all of which function as components of plasma lipoproteins. ApoL-VI (apolipoprotein L-VI), also known as APOE6, is a 343 amino acid protein that localizes to the cytoplasm and belongs to the apolipoprotein L family. Expressed in a variety of tissues, including liver, heart, uterus, spleen, colon, spinal cord, placenta, prostate and mammary gland, apoL-VI is thought to affect the movement of lipids in the cytoplasm and may allow the binding of lipids to organelles. Overproduction of apoL-VI induces apoptosis, suggesting that apoL-VI may also be involved in tumor progression. Like other members of the apolipoprotein L family, apoL-VI is thought to be involved in the development of schizophrenia.

REFERENCES

1. Dunham, I., Shimizu, N., Roe, B.A., Chisoe, S., Hunt, A.R., Collins, J.E., Bruskewich, R., Beare, D.M., Clamp, M., Smink, L.J., Ainscough, R., Almeida, J.P., Babbage, A., Bagguley, C., Bailey, J., Barlow, K., et al. 1999. The DNA sequence of human chromosome 22. *Nature* 402: 489-495.
2. Page, N.M., Butlin, D.J., Lomthaisong, K. and Lowry, P.J. 2001. The human apolipoprotein L gene cluster: identification, classification, and sites of distribution. *Genomics* 74: 71-78.
3. Monajemi, H., Fontijn, R.D., Pannekoek, H. and Horrevoets, A.J. 2002. The apolipoprotein L gene cluster has emerged recently in evolution and is expressed in human vascular tissue. *Genomics* 79: 539-546.
4. Mimmack, M.L., Ryan, M., Baba, H., Navarro-Ruiz, J., Iritani, S., Faull, R.L., McKenna, P.J., Jones, P.B., Arai, H., Starkey, M., Emson, P.C. and Bahn, S. 2002. Gene expression analysis in schizophrenia: reproducible up-regulation of several members of the apolipoprotein L family located in a high-susceptibility locus for schizophrenia on chromosome 22. *Proc. Natl. Acad. Sci. USA* 99: 4680-4685.
5. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607256. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Liu, Z., Lu, H., Jiang, Z., Pastuszyn, A. and Hu, C.A. 2005. Apolipoprotein I6, a novel proapoptotic Bcl-2 homology 3-only protein, induces mitochondria-mediated apoptosis in cancer cells. *Mol. Cancer Res.* 3: 21-31.
7. Liu, Y.L., Fann, C.S., Liu, C.M., Chen, W.J., Wu, J.Y., Hung, S.I., Chen, C.H., Jou, Y.S., Liu, S.K., Hwang, T.J., Hsieh, M.H., Chang, C.C., Yang, W.C., Lin, J.J., Chou, F.H., Faraone, S.V., Tsuang, M.T. and Hwu, H.G. 2008. RASD2, MYH9, and CACNG2 genes at chromosome 22q12 associated with the subgroup of schizophrenia with non-deficit in sustained attention and executive function. *Biol. Psychiatry* 64: 789-796.

CHROMOSOMAL LOCATION

Genetic locus: ApoL6 (mouse) mapping to 15 D3.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

apoL-VI siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see apoL-VI shRNA Plasmid (m): sc-141162-SH and apoL-VI shRNA (m) Lentiviral Particles: sc-141162-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

apoL-VI siRNA (m) is recommended for the inhibition of apoL-VI expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor apoL-VI gene expression knockdown using RT-PCR Primer: apoL-VI (m)-PR: sc-141162-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.