

# ARHGEF19 siRNA (m): sc-141227

## BACKGROUND

Rho GTPases, which play fundamental roles in numerous cellular processes, are initiated by external stimuli that signal through G protein-coupled receptors. ARHGEF19 (Rho guanine nucleotide exchange factor (GEF) 19), also known as WGEF (weakly similar to Rho GEF), is an 802 amino acid protein that contains one DH (DBL-homology) domain, one PH domain and one SH3 domain. Conserved in chimpanzee, canine, bovine, mouse, rat and zebrafish, ARHGEF19 exists as two alternatively spliced isoforms and shares 48% similarity with ARHGEF5. Highly expressed in intestine, with lower expression in liver, heart and kidney, ARHGEF19 displays GEF activity towards Rho A family members, with possible activity against Rac 1 and Cdc42. ARHGEF19 regulates adipogenesis by modifying DNA methylation. Overexpression of ARHGEF19 may be linked to filopodia and stress fiber formation, membrane ruffling and lamellipodia.

## REFERENCES

- Wang, Y., Suzuki, H., Yokoo, T., Tada-Iida, K., Kihara, R., Miura, M., Watanabe, K., Sone, H., Shimano, H., Toyoshima, H. and Yamada, N. 2004. WGEF is a novel Rho GEF expressed in intestine, liver, heart, and kidney. *Biochem. Biophys. Res. Commun.* 324: 1053-1058.
- Sahin, M., Greer, P.L., Lin, M.Z., Poucher, H., Eberhart, J., Schmidt, S., Wright, T.M., Shamah, S.M., O'Connell, S., Cowan, C.W., Hu, L., Goldberg, J.L., Debant, A., Corfas, G., Krull, C.E. and Greenberg, M.E. 2005. Eph-dependent tyrosine phosphorylation of ephexin1 modulates growth cone collapse. *Neuron* 46: 191-204.
- Nowicki, M., Kosacka, J., Brossmer, R., Spaniel-Borowski, K. and Borlak, J. 2007. The myelin-associated glycoprotein inhibitor BENZ induces outgrowth and survival of rat dorsal root ganglion cell cultures. *J. Neurosci. Res.* 85: 3053-3063.
- McKown, R.L., Wang, N., Raab, R.W., Karnati, R., Zhang, Y., Williams, P.B. and Laurie, G.W. 2009. Lacritin and other new proteins of the lacrimal functional unit. *Exp. Eye Res.* 88: 848-858.
- Mikkelsen, T., Brodie, C., Finniss, S., Berens, M.E., Rennert, J.L., Nelson, K., Lemke, N., Brown, S.L., Hahn, D., Neuteboom, B. and Goodman, S.L. 2009. Radiation sensitization of glioblastoma by cilengitide has unanticipated schedule-dependency. *Int. J. Cancer* 124: 2719-2727.
- Horii, T., Morita, S., Kimura, M. and Hatada, I. 2009. Epigenetic regulation of adipocyte differentiation by a Rho guanine nucleotide exchange factor, WGEF. *PLoS ONE* 4: e5809.
- Cordeaux, Y., Tattersall, M., Charnock-Jones, D.S. and Smith, G.C. 2010. Effects of medroxyprogesterone acetate on gene expression in myometrial explants from pregnant women. *J. Clin. Endocrinol. Metab.* 95: E437-E447.
- Hiramoto-Yamaki, N., Takeuchi, S., Ueda, S., Harada, K., Fujimoto, S., Negishi, M. and Katoh, H. 2010. Ephexin4 and EphA2 mediate cell migration through a RhoG-dependent mechanism. *J. Cell Biol.* 190: 461-477.
- Loirand, G. and Pacaud, P. 2010. The role of Rho protein signaling in hypertension. *Nat. Rev. Cardiol.* 7: 637-647.

## CHROMOSOMAL LOCATION

Genetic locus: Arhgef19 (mouse) mapping to 4 D3.

## PRODUCT

ARHGEF19 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ARHGEF19 shRNA Plasmid (m): sc-141227-SH and ARHGEF19 shRNA (m) Lentiviral Particles: sc-141227-V as alternate gene silencing products.

For independent verification of ARHGEF19 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-141227A and sc-141227B.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

ARHGEF19 siRNA (m) is recommended for the inhibition of ARHGEF19 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ARHGEF19 gene expression knockdown using RT-PCR Primer: ARHGEF19 (m)-PR: sc-141227-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.