

ASPHD1 siRNA (m): sc-141305

BACKGROUND

ASPHD1 (aspartate beta-hydroxylase domain containing 1) is a 390 amino acid single-pass type II membrane protein belonging to the aspartyl/asparaginyl β -hydroxylase family. Conserved in chimpanzee, canine, bovine, rat and zebrafish, ASPHD1 participates in oxidoreductase activity, acting on single donors with incorporation of two atoms of oxygen, and peptide-aspartate β -dioxygenase activity. Encoded by a gene that maps to human chromosome 16p11.2, ASPHD1 is one of 25 genes in a microdeletion associated with endocardial fibroelastosis. Duplication in 16p11.2 is linked to malignant migrating partial seizure disorder of infancy, and recurrent rearrangements in 16p11.2 are linked to numerous neurodevelopmental disorders, including autism, mental retardation and schizophrenia. Paroxysmal kinesigenic choreoathetosis (PKC) is also linked to 16p11-q21 mutations.

REFERENCES

- Gerhard, D.S., Wagner, L., Feingold, E.A., Shenmen, C.M., Grouse, L.H., Schuler, G., Klein, S.L., Old, S., Rasooly, R., Good, P., Guyer, M., Peck, A.M., Derge, J.G., Lipman, D., Collins, F.S., Jang, W., Sherry, S., et al. 2004. The status, quality, and expansion of the NIH full-length cDNA project: the Mammalian Gene Collection (MGC). *Genome Res.* 14: 2121-2127.
- Kikuchi, T., Nomura, M., Tomita, H., Harada, N., Kanai, K., Konishi, T., Yasuda, A., Matsuura, M., Kato, N., Yoshiura, K. and Niikawa, N. 2007. Paroxysmal kinesigenic choreoathetosis (PKC): confirmation of linkage to 16p11-q21, but unsuccessful detection of mutations among 157 genes at the PKC-critical region in seven PKC families. *J. Hum. Genet.* 52: 334-341.
- Pollard, P.J., Loenarz, C., Mole, D.R., McDonough, M.A., Gleadle, J.M., Schofield, C.J. and Ratcliffe, P.J. 2008. Regulation of Jumonji-domain-containing histone demethylases by hypoxia-inducible factor (HIF)-1 α . *Biochem. J.* 416: 387-394.
- Daub, H., Olsen, J.V., Bairlein, M., Gnad, F., Oppermann, F.S., Körner, R., Greff, Z., Keri, G., Stemmann, O. and Mann, M. 2008. Kinase-selective enrichment enables quantitative phosphoproteomics of the kinome across the cell cycle. *Mol. Cell* 31: 438-448.
- Bedoyan, J.K., Kumar, R.A., Sudi, J., Silverstein, F., Ackley, T., Iyer, R.K., Christian, S.L. and Martin, D.M. 2010. Duplication 16p11.2 in a child with infantile seizure disorder. *Am. J. Med. Genet. A* 152A: 1567-1574.
- Puvabanditsin, S., Nagar, M.S., Joshi, M., Lambert, G., Garrow, E. and Brandsma, E. 2010. Microdeletion of 16p11.2 associated with endocardial fibroelastosis. *Am. J. Med. Genet. A* 152A: 2383-2386.
- Hsu, P.Y., Hsu, H.K., Singer, G.A., Yan, P.S., Rodriguez, B.A., Liu, J.C., Weng, Y.I., Deatherage, D.E., Chen, Z., Pereira, J.S., Lopez, R., Russo, J., Wang, Q., Lamartiniere, C.A., Nephew, K.P. and Huang, T.H. 2010. Estrogen-mediated epigenetic repression of large chromosomal regions through DNA looping. *Genome Res.* 20: 733-744.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: *Asphd1* (mouse) mapping to 7 F3.

PRODUCT

ASPHD1 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ASPHD1 shRNA Plasmid (m): sc-141305-SH and ASPHD1 shRNA (m) Lentiviral Particles: sc-141305-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ASPHD1 siRNA (m) is recommended for the inhibition of ASPHD1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ASPHD1 gene expression knockdown using RT-PCR Primer: ASPHD1 (m)-PR: sc-141305-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.