

BC018242 siRNA (m): sc-141514

BACKGROUND

Phosphatidate phosphatases are a family of integral membrane glycoproteins that dephosphorylate a variety of lipid phosphates and play a role in signal transduction via the phospholipase D pathway. PAP-2 proteins function independently of Mg^{2+} and are insensitive to NEM (N-ethylmaleimide) inhibition. The lipid phosphates degraded by this family include ceramide 1-phosphate (C1P), sphingosine 1-phosphate (S1P), phosphatidic acid (PA) and lysophosphatidic acid (LPA). LPPR2 (lipid phosphate phosphatase-related protein type 2), also known as PRG4 (plasticity-related gene 4 protein), is a 343 amino acid multi-pass membrane protein that mediates lipid phosphate phosphatase activity and may be involved in mechanisms of neuronal plasticity. There are two isoforms of LPPR2 that are produced as a result of alternative splicing events.

REFERENCES

- Bräuer, A.U., Savaskan, N.E., Kühn, H., Pohn, S., Ninnemann, O. and Nitsch, R. 2003. A new phospholipid phosphatase, PRG-1, is involved in axon growth and regenerative sprouting. *Nat. Neurosci.* 6: 572-578.
- Savaskan, N.E., Bräuer, A.U. and Nitsch, R. 2004. Molecular cloning and expression regulation of PRG-3, a new member of the plasticity-related gene family. *Eur. J. Neurosci.* 19: 212-220.
- Sigal, Y.J., McDermott, M.I. and Morris, A.J. 2005. Integral membrane lipid phosphatases/phosphotransferases: common structure and diverse functions. *Biochem. J.* 387: 281-293.
- Tanic, N., Brkic, G., Dimitrijevic, B., Dedovic-Tanic, N., Gefen, N., Benharroch, D. and Gopas, J. 2006. Identification of differentially expressed mRNA transcripts in drug-resistant versus parental human melanoma cell lines. *Anticancer Res.* 26: 2137-2142.
- Theofilopoulos, S., Lykidis, A., Leondaritis, G. and Mangoura, D. 2008. Novel function of the human presqualene diphosphate phosphatase as a type II phosphatidate phosphatase in phosphatidylcholine and triacylglyceride biosynthesis pathways. *Biochim. Biophys. Acta* 1781: 731-742.
- Brindley, D.N., Pilquil, C., Sariahmetoglu, M. and Reue, K. 2009. Phosphatidate degradation: phosphatidate phosphatases (lipins) and lipid phosphate phosphatases. *Biochim. Biophys. Acta* 1791: 956-961.
- Brindley, D.N. and Pilquil, C. 2009. Lipid phosphate phosphatases and signaling. *J. Lipid Res.* 50: S225-S230.
- Broggini, T., Nitsch, R. and Savaskan, N.E. 2010. Plasticity-related gene 5 (PRG5) induces filopodia and neurite growth and impedes lysophosphatidic acid- and nogo-A-mediated axonal retraction. *Mol. Biol. Cell* 21: 521-537.

CHROMOSOMAL LOCATION

Genetic locus: BC018242 (mouse) mapping to 9 A3.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

BC018242 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see BC018242 shRNA Plasmid (m): sc-141514-SH and BC018242 shRNA (m) Lentiviral Particles: sc-141514-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

BC018242 siRNA (m) is recommended for the inhibition of BC018242 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor BC018242 gene expression knockdown using RT-PCR Primer: BC018242 (m)-PR: sc-141514-PR (20 μ l). Annealing temperature for the primers should be $55-60^{\circ}$ C and the extension temperature should be $68-72^{\circ}$ C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.