

CEP152 siRNA (m): sc-142278

BACKGROUND

CEP152, also known as centrosomal protein of 152 kDa, MCPH4 or SCKL5, is a 1,654 amino acid protein that is required for centrosome duplication. CEP152 is responsible for regulating genomic integrity as well as initiating ATR-mediated checkpoint signaling in response to DNA damage. CEP152 interacts with Sak and CENPJ, which are important for centriole formation as well as the cyclin-dependent kinase 2-interacting protein (CINP). Defects in CEP152 leads to Seckel syndrome type 5 (SCKL5) and microcephaly primary type 4 (MCPH4), both of which are congenital diseases and cause mental retardation. CEP152 exists as three different alternatively spliced isoforms and is post-translationally phosphorylated at amino acid residue 1,258 (tyrosine) and 1,405 (serine). The gene encoding CEP152 maps to human chromosome 15.

REFERENCES

1. Nagase, T., Ishikawa, K., Suyama, M., Kikuno, R., Hirosawa, M., Miyajima, N., Tanaka, A., Kotani, H., Nomura, N. and Ohara, O. 1998. Prediction of the coding sequences of unidentified human genes. XII. The complete sequences of 100 new cDNA clones from brain which code for large proteins *in vitro*. DNA Res. 5: 355-364.
2. Andersen, J.S., Wilkinson, C.J., Mayor, T., Mortensen, P., Nigg, E.A. and Mann, M. 2003. Proteomic characterization of the human centrosome by protein correlation profiling. Nature 426: 570-574.
3. Dephoure, N., Zhou, C., Villen, J., Beausoleil, S.A., Bakalarski, C.E., Elledge, S.J. and Gygi, S.P. 2008. A quantitative atlas of mitotic phosphorylation. Proc. Natl. Acad. Sci. USA 105: 10762-10767.
4. Heibeck, T.H., Ding, S.J., Opresko, L.K., Zhao, R., Schepmoes, A.A., Yang, F., Tolmachev, A.V., Monroe, M.E., Camp, D.G., Smith, R.D., Wiley, H.S. and Qian, W.J. 2009. An extensive survey of tyrosine phosphorylation revealing new sites in human mammary epithelial cells. J. Proteome Res. 8: 3852-3861.
5. Hatch, E.M., Kulukian, A., Holland, A.J., Cleveland, D.W. and Stearns, T. 2010. Cep152 interacts with Plk4 and is required for centriole duplication. J. Cell Biol. 191: 721-729.
6. Cizmecioglu, O., Arnold, M., Bahtz, R., Settele, F., Ehret, L., Haselmann-Weiss, U., Antony, C. and Hoffmann, I. 2010. Cep152 acts as a scaffold for recruitment of Plk4 and CPAP to the centrosome. J. Cell Biol. 191: 731-739.
7. Dzhindzhev, N.S., Yu, Q.D., Weiskopf, K., Tzolovsky, G., Cunha-Ferreira, I., Riparbelli, M., Rodrigues-Martins, A., Bettencourt-Dias, M., Callaini, G. and Glover, D.M. 2010. Asterless is a scaffold for the onset of centriole assembly. Nature 467: 714-718.
8. Kalay, E., Yigit, G., Aslan, Y., Brown, K.E., Pohl, E., Bicknell, L.S., Kayserili, H., Li, Y., Tüysüz, B., Nürnberg, G., Kiess, W., Koegl, M., Baessmann, I., Buruk, K., Toraman, B., Kayipmaz, S., Kul, S., Ikbali, M., Turner, D.J., et al. 2011. CEP152 is a genome maintenance protein disrupted in Seckel syndrome. Nat. Genet. 43: 23-26.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: Cep152 (mouse) mapping to 2 F1.

PRODUCT

CEP152 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CEP152 shRNA Plasmid (m): sc-142278-SH and CEP152 shRNA (m) Lentiviral Particles: sc-142278-V as alternate gene silencing products.

For independent verification of CEP152 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-142278A, sc-142278B and sc-142278C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

CEP152 siRNA (m) is recommended for the inhibition of CEP152 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CEP152 gene expression knockdown using RT-PCR Primer: CEP152 (m)-PR: sc-142278-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.