

Clpb siRNA (m): sc-142401

BACKGROUND

HSP 78, also known as CLPB (caseinolytic peptidase B protein homolog) or SKD3 (suppressor of potassium transport defect 3), is a 707 amino acid protein that belongs to the clpA/clpB family. Containing four ANK repeats, HSP 78 may be involved in the regulation of ATPase and the secretion/protein trafficking process. Existing as two alternatively spliced isoforms, HSP 78 maps to human chromosome 11q13.4 and mouse chromosome 7 E3. Human chromosome 11 houses over 1,400 genes, comprises nearly 4% of the human genome and is considered a gene and disease association dense chromosome. Jervell and Lange-Nielsen syndrome, Jacobsen syndrome, Niemann-Pick disease, hereditary angioedema and Smith-Lemli-Opitz syndrome are associated with defects in genes that maps to chromosome 11.

REFERENCES

1. Perier, F., Radeke, C.M., Raab-Graham, K.F. and Vandenberg, C.A. 1995. Expression of a putative ATPase suppresses the growth defect of a yeast potassium transport mutant: identification of a mammalian member of the Clp/HSP104 family. *Gene* 152: 157-163.
2. Rottgers, K., Zufall, N., Guiard, B. and Voos, W. 2002. The ClpB homolog Hsp78 is required for the efficient degradation of proteins in the mitochondrial matrix. *J. Biol. Chem.* 277: 45829-45837.
3. Jira, P.E., Waterham, H.R., Wanders, R.J., Smeitink, J.A., Sengers, R.C. and Wevers, R.A. 2003. Smith-Lemli-Opitz syndrome and the DHCR7 gene. *Ann. Hum. Genet.* 67: 269-280.
4. Lewandowska, A., Gierszewska, M., Marszalek, J. and Liberek, K. 2006. Hsp78 chaperone functions in restoration of mitochondrial network following heat stress. *Biochim. Biophys. Acta* 1763: 141-151.
5. Bösl, B., Grimminger, V. and Walter, S. 2006. The molecular chaperone Hsp104—a molecular machine for protein disaggregation. *J. Struct. Biol.* 156: 139-148.
6. Leidhold, C., von Janowsky, B., Becker, D., Bender, T. and Voos, W. 2006. Structure and function of Hsp78, the mitochondrial ClpB homolog. *J. Struct. Biol.* 156: 149-164.
7. Schuchman, E.H. 2007. The pathogenesis and treatment of acid sphingomyelinase-deficient Niemann-Pick disease. *J. Inherit. Metab. Dis.* 30: 654-663.
8. Bhuiyan, Z.A., Momenah, T.S., Amin, A.S., Al-Khadra, A.S., Alders, M., Wilde, A.A. and Mannens, M.M. 2008. An intronic mutation leading to incomplete skipping of exon-2 in KCNQ1 rescues hearing in Jervell and Lange-Nielsen syndrome. *Prog. Biophys. Mol. Biol.* 98: 319-327.
9. Desantis, M.E. and Shorter, J. 2012. The elusive middle domain of Hsp104 and ClpB: Location and function. *Biochim Biophys Acta* 1823: 29-39.

CHROMOSOMAL LOCATION

Genetic locus: Clpb (mouse) mapping to 7 E3.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Clpb siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Clpb shRNA Plasmid (m): sc-142401-SH and Clpb shRNA (m) Lentiviral Particles: sc-142401-V as alternate gene silencing products.

For independent verification of Clpb (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-142401A, sc-142401B and sc-142401C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Clpb siRNA (m) is recommended for the inhibition of Clpb expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Clpb gene expression knockdown using RT-PCR Primer: Clpb (m)-PR: sc-142401-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.