

# Cml3 siRNA (m): sc-142414

## BACKGROUND

Cml3 (camello-like protein 3), also known as probable N-acetyltransferase CML3, is a 226 amino acid multi-pass membrane protein that contains one N-acetyltransferase domain and belongs to the camello family. Cml3 may play a role in regulation of gastrulation, and is encoded by a gene that maps to mouse chromosome 6 C3. The human homolog of Cml3, known as NAT-8L, N-acetyltransferase 8-like protein, N-acetylaspartate synthetase or CML3 (camello-like protein 3), is a 302 amino acid protein that plays a role in the regulation of lipogenesis by producing N-acetylaspartate acid. Defects in NAT-8L are the cause of N-acetylaspartate deficiency (NACED), a metabolic disorder resulting in truncal ataxia, marked developmental delay, seizures and secondary microcephaly.

## REFERENCES

1. Popsueva, A.E., Luchinskaya, N.N., Ludwig, A.V., Zinovjeva, O.Y., Poteryaev, D.A., Feigelman, M.M., Ponomarev, M.B., Berekelya, L. and Belyavsky, A.V. 2001. Overexpression of camello, a member of a novel protein family, reduces blastomere adhesion and inhibits gastrulation in *Xenopus laevis*. *Dev. Biol.* 234: 483-496.
2. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 610647. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Arun, P., Moffett, J.R. and Namboodiri, A.M. 2009. Evidence for mitochondrial and cytoplasmic N-acetylaspartate synthesis in SH-SY5Y neuroblastoma cells. *Neurochem. Int.* 55: 219-225.
4. Wiame, E., Tyteca, D., Pierrot, N., Collard, F., Amyere, M., Noel, G., Desmedt, J., Nassogne, M.C., Vikkula, M., Octave, J.N., Vincent, M.F., Courtoy, P.J., Boltshauser, E. and van Schaftingen, E. 2010. Molecular identification of aspartate N-acetyltransferase and its mutation in hypoaetylaspartia. *Biochem. J.* 425: 127-136.
5. Ariyannur, P.S., Moffett, J.R., Manickam, P., Pattabiraman, N., Arun, P., Nitta, A., Nabeshima, T., Madhavarao, C.N. and Namboodiri, A.M. 2010. Methamphetamine-induced neuronal protein NAT8L is the NAA biosynthetic enzyme: implications for specialized acetyl coenzyme A metabolism in the CNS. *Brain Res.* 1335: 1-13.
6. Becker, I., Lodder, J., Gieselmann, V. and Eckhardt, M. 2010. Molecular characterization of N-acetylasparylglutamate synthetase. *J. Biol. Chem.* 285: 29156-29164.
7. Tahay, G., Wiame, E., Tyteca, D., Courtoy, P.J. and Van Schaftingen, E. 2012. Determinants of the enzymatic activity and the subcellular localization of aspartate N-acetyltransferase. *Biochem. J.* 441: 105-112.

## CHROMOSOMAL LOCATION

Genetic locus: Cml3 (mouse) mapping to 6 C3.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

Cml3 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Cml3 shRNA Plasmid (m): sc-142414-SH and Cml3 shRNA (m) Lentiviral Particles: sc-142414-V as alternate gene silencing products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

Cml3 siRNA (m) is recommended for the inhibition of Cml3 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Cml3 gene expression knockdown using RT-PCR Primer: Cml3 (m)-PR: sc-142414-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.