



copine 5 shRNA (m) Lentiviral Particles: sc-142506-V

BACKGROUND

Copine 5, also known as CPNE5 or CPN5, is a 593 amino acid calcium-dependent membrane-binding protein that contains two N-terminal C2 domains and one WWFA domain. Expressed as multiple alternatively spliced isoforms, copine 5 is thought to play a role in membrane trafficking, specifically exhibiting calcium-dependent phospholipid binding properties. The gene encoding copine 5 maps to human chromosome 6, which contains 170 million base pairs and comprises nearly 6% of the human genome. Deletion of a portion of the q arm of chromosome 6 is associated with early onset intestinal cancer, suggesting the presence of a cancer susceptibility locus. Additionally, Porphyria cutanea tarda, Parkinson's disease, Stickler syndrome and a susceptibility to bipolar disorder are all associated with genes that map to chromosome 6.

REFERENCES

1. Tripodis, N., Mason, R., Humphray, S.J., Davies, A.F., Herberg, J.A., Trowsdale, J., Nizetic, D., Senger, G. and Ragoussis, J. 1998. Physical map of human 6p21.2-6p21.3: region flanking the centromeric end of the major histocompatibility complex. *Genome Res.* 8: 631-643.
2. Creutz, C.E., Tomsig, J.L., Snyder, S.L., Gautier, M.C., Skouri, F., Beisson, J. and Cohen, J. 1998. The copines, a novel class of C2 domain-containing, calcium-dependent, phospholipid-binding proteins conserved from *Paramecium* to humans. *J. Biol. Chem.* 273: 1393-1402.
3. Online Mendelian Inheritance in Man, OMIM™. 1999. Johns Hopkins University, Baltimore, MD. MIM Number: 604209. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Vuoristo, M.M., Pappas, J.G., Jansen, V. and Ala-Kokko, L. 2004. A stop codon mutation in COL11A2 induces exon skipping and leads to non-ocular Stickler syndrome. *Am. J. Med. Genet. A* 130A: 160-164.

CHROMOSOMAL LOCATION

Genetic locus: Cpne5 (mouse) mapping to 17 A3.3.

PRODUCT

copine 5 shRNA (m) Lentiviral Particles are concentrated, transduction-ready viral particles containing a target-specific construct that encodes a 19-25 nt (plus hairpin) shRNA designed to knock down gene expression. Each vial contains 200 µl frozen stock containing 1.0×10^6 infectious units of virus (IFU) in Dulbecco's Modified Eagle's Medium with 25 mM HEPES pH 7.3. Suitable for 10-20 transductions. Also see copine 5 siRNA (m): sc-142506 and copine 5 shRNA Plasmid (m): sc-142506-SH as alternate gene silencing products.

APPLICATIONS

copine 5 shRNA (m) Lentiviral Particles is recommended for the inhibition of copine 5 expression in mouse cells.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

SUPPORT REAGENTS

Control shRNA Lentiviral Particles: sc-108080. Available as 200 µl frozen viral stock containing 1.0×10^6 infectious units of virus (IFU); contains an shRNA construct encoding a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor copine 5 gene expression knockdown using RT-PCR Primer: copine 5 (m)-PR: sc-142506-PR (20 µl). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

BIOSAFETY

Lentiviral particles can be employed in standard Biosafety Level 2 tissue culture facilities (and should be treated with the same level of caution as with any other potentially infectious reagent). Lentiviral particles are replication-incompetent and are designed to self-inactivate after transduction and integration of shRNA constructs into genomic DNA of target cells.

STORAGE

Store lentiviral particles at -80° C. Stable for at least one year from the date of shipment. Once thawed, particles can be stored at 4° C for up to one week. Avoid repeated freeze thaw cycles.

RESEARCH USE

The purchase of this product conveys to the buyer the nontransferable right to use the purchased amount of the product and all replicates and derivatives for research purposes conducted by the buyer in his laboratory only (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party, or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes.