

COX8c siRNA (m): sc-142535

BACKGROUND

The cytochrome c oxidase (COX) family of proteins function as the final electron donor in the respiratory chain to drive a proton gradient across the inner mitochondrial membrane, ultimately resulting in the production of water and ATP. The mammalian COX apoenzyme is a dimer, with each monomer consisting of 13 subunits, some of which are mitochondrial and some of which are nuclear. The COX8 (cytochrome c oxidase subunit VIII) subunits are nuclear and have muscle and non-muscle-specific isoforms. COX8 exists as three isoforms: COX8a, a liver and heart isoform, Cox8b, a heart-specific isoform, and Cox8c, whose expression pattern has yet to be elucidated. All three Cox8 isoforms exist as components of the COX complex and play an important role in electron transport.

REFERENCES

1. Patterson, T.E. and Poyton, R.O. 1986. COX8, the structural gene for yeast cytochrome c oxidase subunit VIII. DNA sequence and gene disruption indicate that subunit VIII is required for maximal levels of cellular respiration and is derived from a precursor which is extended at both its NH₂ and COOH termini. *J. Biol. Chem.* 261: 17192-17197.
2. Rizzuto, R., Nakase, H., Darras, B., Francke, U., Fabrizi, G.M., Mengel, T., Walsh, F., Kadenbach, B., DiMauro, S. and Schon, E.A. 1989. A gene specifying subunit VIII of human cytochrome c oxidase is localized to chromosome 11 and is expressed in both muscle and non-muscle tissues. *J. Biol. Chem.* 264: 10595-10600.
3. Bonne, G., Carrier, L., Schwartz, K. and Komajda, M. 1995. The COX8 gene is not the disease gene of the CMH4 locus in familial hypertrophic cardiomyopathy. *J. Med. Genet.* 32: 670-671.
4. Lomax, M.I., Riggs, P.K. and Womack, J.E. 1995. Structure and chromosomal location of the bovine gene for the heart muscle isoform of cytochrome c oxidase subunit VIII. *Mamm. Genome* 6: 118-122.
5. Hüttemann, M., Schmidt, T.R. and Grossman, L.I. 2003. A third isoform of cytochrome c oxidase subunit VIII is present in mammals. *Gene* 312: 95-102.
6. Khalimonchuk, O. and Rödel, G. 2005. Biogenesis of cytochrome c oxidase. *Mitochondrion* 5: 363-388.
7. Fontanesi, F., Soto, I.C. and Barrientos, A. 2008. Cytochrome c oxidase biogenesis: new levels of regulation. *IUBMB Life* 60: 557-568.
8. Barrientos, A., Gouget, K., Horn, D., Soto, I.C. and Fontanesi, F. 2009. Suppression mechanisms of COX assembly defects in yeast and human: insights into the COX assembly process. *Biochim. Biophys. Acta* 1793: 97-107.

CHROMOSOMAL LOCATION

Genetic locus: Cox8c (mouse) mapping to 12 E.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

COX8c siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see COX8c shRNA Plasmid (m): sc-142535-SH and COX8c shRNA (m) Lentiviral Particles: sc-142535-V as alternate gene silencing products.

For independent verification of COX8c (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-142535A and sc-142535B.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

COX8c siRNA (m) is recommended for the inhibition of COX8c expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RESEARCH USE

For research use only, not for use in diagnostic procedures.