

# CPA4 siRNA (m): sc-142538

## BACKGROUND

CPA4 (carboxypeptidase A4), also known as CPA3 or (carboxypeptidase A3), is a 421 amino acid protein that functions as a metalloprotease and is likely involved in the histone hyperacetylation pathway. A secreted protein and member of the peptidase M14 family, CPA4 is expressed in fetal lung, brain, heart, adrenal gland, intestine, kidney and liver. CPA4 was first identified in PC-3, a prostate cancer cell line, and is also expressed in adult benign hyper-trophic prostate. CPA4 is considered a potential marker for aggressive prostate cancer and is upregulated by inhibitors of histone deacetylation. CPA4 binds one zinc ion per subunit and is encoded by a gene that maps to human chromosome 7q32.2. Chromosome 7 houses over 1,000 genes, comprises nearly 5% of the human genome and has been linked to Osteogenesis imperfecta, Pendred syndrome, Lissencephaly, Citrullinemia and Shwachman-Diamond syndrome.

## REFERENCES

1. Tsiouras, P., Myers, J.C., Ramirez, F. and Prockop, D.J. 1983. Restriction  $\alpha$ human type I procollagen. Application to a family with an autosomal dominant form of osteogenesis imperfecta. *J. Clin. Invest.* 72: 1262-1267.
2. Huang, H., Reed, C.P., Zhang, J.S., Shridhar, V., Wang, L. and Smith, D.I. 1999. Carboxypeptidase A3 (CPA3): a novel gene highly induced by histone deacetylase inhibitors during differentiation of prostate epithelial cancer cells. *Cancer Res.* 59: 2981-2988.
3. Witte, J.S., Goddard, K.A., Conti, D.V., Elston, R.C., Lin, J., Suarez, B.K., Broman, K.W., Burmester, J.K., Weber, J.L. and Catalona, W.J. 2000. Genomewide scan for prostate cancer-aggressiveness loci. *Am. J. Hum. Genet.* 67: 92-99.
4. Iwasaki, S., Usami, S., Abe, S., Isoda, H., Watanabe, T. and Hoshino, T. 2001. Long-term audiological feature in Pendred syndrome caused by PDS mutation. *Arch. Otolaryngol. Head Neck Surg.* 127: 705-708.
5. Kayashima, T., Yamasaki, K., Yamada, T., Sakai, H., Miwa, N., Ohta, T., Yoshiura, K., Matsumoto, N., Nakane, Y., Kanetake, H., Ishino, F., Niiikawa, N. and Kishino, T. 2003. The novel imprinted carboxypeptidase A4 gene (CPA4) in the 7q32 imprinting domain. *Hum. Genet.* 112: 220-226.
6. Bentley, L., Nakabayashi, K., Monk, D., Beechey, C., Peters, J., Birjandi, Z., Khayat, F.E., Patel, M., Preece, M.A., Stanier, P., Scherer, S.W. and Moore, G.E. 2003. The imprinted region on human chromosome 7q32 extends to the carboxypeptidase A gene cluster: an imprinted candidate for Silver-Russell syndrome. *J. Med. Genet.* 40: 249-256.
7. Pallarès, I., Bonet, R., García-Castellanos, R., Ventura, S., Aviles, F.X., Vendrell, J. and Gomis-Rüth, F.X. 2005. Structure of human carboxypeptidase A4 with its endogenous protein inhibitor, latexin. *Proc. Natl. Acad. Sci. USA* 102: 3978-3983.
8. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2005. Johns Hopkins University, Baltimore, MD. MIM Number: 607635. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

## CHROMOSOMAL LOCATION

Genetic locus: Cpa4 (mouse) mapping to 6 A3.3.

## PRODUCT

CPA4 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CPA4 shRNA Plasmid (m): sc-142538-SH and CPA4 shRNA (m) Lentiviral Particles: sc-142538-V as alternate gene silencing products.

For independent verification of CPA4 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-142538A, sc-142538B and sc-142538C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

CPA4 siRNA (m) is recommended for the inhibition of CPA4 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CPA4 gene expression knockdown using RT-PCR Primer: CPA4 (m)-PR: sc-142538-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.