



CPVL siRNA (m): sc-142551

BACKGROUND

Carboxypeptidases function as proteases and cleave single amino acids from the C-terminus of peptides or proteins. There are three main groups of carboxypeptidases, namely serine-, cysteine- and metallo-enzymes. CPVL (carboxypeptidase, vitellogenic-like), also known as HVLP (VCP-like protein), is a serine carboxypeptidase that is similar to the vitellogenic carboxypeptidase found in mosquito ovaries. Belonging to the peptidase S10 family, CPVL is expressed in myeloid cells of the immune system and is also found in spleen, kidneys, placenta and heart. CPVL contains four putative N-glycosylation sites and a serine carboxypeptidase active site. During monocyte maturation into macrophages, CPVL expression is induced. This suggests a possible role for CPVL in phagocytosis, antigen processing and organization of the innate immune response.

REFERENCES

1. Mahoney, J.A., Ntoli, B., da Silva, R.P., Gordon, S. and McKnight, A.J. 2001. Cloning and characterization of CPVL, a novel serine carboxypeptidase, from human macrophages. *Genomics* 72: 243-251.
2. Stanton, L.A., Fenhalls, G., Lucas, A., Gough, P., Greaves, D.R., Mahoney, J.A., Helden, P. and Gordon, S. 2003. Immunophenotyping of macrophages in human pulmonary tuberculosis and sarcoidosis. *Int. J. Exp. Pathol.* 84: 289-304.
3. Sleat, D.E., Wang, Y., Sohar, I., Lackland, H., Li, Y., Li, H., Zheng, H. and Lobel, P. 2006. Identification and validation of mannose 6-phosphate glycoproteins in human plasma reveal a wide range of lysosomal and non-lysosomal proteins. *Mol. Cell. Proteomics* 5: 1942-1956.
4. Lee, T.H., Streb, J.W., Georger, M.A. and Miano, J.M. 2006. Tissue expression of the novel serine carboxypeptidase Scep1. *J. Histochem. Cytochem.* 54: 701-711.
5. Mittapalli, O., Wise, I.L. and Shukle, R.H. 2006. Characterization of a serine carboxypeptidase in the salivary glands and fat body of the orange wheat blossom midge, *Sitotiplosis mosellana* (Diptera: Cecidomyiidae). *Insect Biochem. Mol. Biol.* 36: 154-160.
6. Harris, J., Schwinn, N., Mahoney, J.A., Lin, H.H., Shaw, M., Howard, C.J., da Silva, R.P. and Gordon, S. 2006. A vitellogenic-like carboxypeptidase expressed by human macrophages is localized in endoplasmic reticulum and membrane ruffles. *Int. J. Exp. Pathol.* 87: 29-39.

CHROMOSOMAL LOCATION

Genetic locus: Cplv (mouse) mapping to 6 B3.

PRODUCT

CPVL siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CPVL shRNA Plasmid (m): sc-142551-SH and CPVL shRNA (m) Lentiviral Particles: sc-142551-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

CPVL siRNA (m) is recommended for the inhibition of CPVL expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CPVL gene expression knockdown using RT-PCR Primer: CPVL (m)-PR: sc-142551-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.