



EF-G1 siRNA (m): sc-143307

BACKGROUND

Protein translation consists of initiation, elongation and termination stages that are catalyzed by several protein factors. In mitochondria, the elongation stage requires three elongation factors: Ts, Tu and G. EF-G1 (elongation factor G1) is a 751 amino acid protein containing a GTP-binding elongation factor signature and a GTP-binding domain. This mitochondrial protein promotes the GTP-dependent translocation of the nascent protein chain from the ribosomal A-site to the P-site. Defects in the gene encoding EF-G1 lead to a mutated GTP-binding domain is the cause of combined oxidative phosphorylation deficiency (COXPD), a disease that results in early fatal progressive hepatoenkephalopathy.

REFERENCES

1. Gao, J., et al. 2001. Cloning and characterization of human and mouse mitochondrial elongation factor G, GFM and Gfm, and mapping of GFM to human chromosome 3q25.1-q26.2. *Genomics* 74: 109-114.
2. Hammarsund, M., et al. 2001. Identification and characterization of two novel human mitochondrial elongation factor genes, hEFG2 and hEFG1, phylogenetically conserved through evolution. *Hum. Genet.* 109: 542-550.
3. Zeviani, M., et al. 2004. Mitochondrial disorders. *Brain* 127: 2153-2172.
4. Coenen, M.J., et al. 2004. Mutant mitochondrial elongation factor G1 and combined oxidative phosphorylation deficiency. *N. Engl. J. Med.* 351: 2080-2086.
5. Smeitink, J.A., et al. 2006. Distinct clinical phenotypes associated with a mutation in the mitochondrial translation elongation factor EFTs. *Am. J. Hum. Genet.* 79: 869-877.
6. Valente, L., et al. 2007. Infantile encephalopathy and defective mitochondrial DNA translation in patients with mutations of mitochondrial elongation factors EFG1 and EFTu. *Am. J. Hum. Genet.* 80: 44-58.
7. Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 606639. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: Gfm1 (mouse) mapping to 3 E1.

PRODUCT

EF-G1 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see EF-G1 shRNA Plasmid (m): sc-143307-SH and EF-G1 shRNA (m) Lentiviral Particles: sc-143307-V as alternate gene silencing products.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

EF-G1 siRNA (m) is recommended for the inhibition of EF-G1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor EF-G1 gene expression knockdown using RT-PCR Primer: EF-G1 (m)-PR: sc-143307-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.