

Elastase-3 siRNA (m): sc-144629

BACKGROUND

Elastases belong to a subfamily of serine proteases that function to catalytically hydrolyze a wide variety of proteins. Elastase-3, also known as Ela3, Ela3b or Cela3b (chymotrypsin-like elastase family, member 3B), is a 255 amino acid murine protein that belongs to the peptidase S1 family and is the homolog of human Elastase-3B. There are six human Elastase genes that encode structurally similar proteins, Elastase-1, 2, 2A, 2B, 3A, and 3B. Unlike the other Elastase proteins, Elastase-3B has little elastolytic activity. Secreted as a zymogen from the pancreas, Elastase-3B possesses similar functions as Elastase-3A, in which it acts as a digestive enzyme in the intestine and participates in metabolic degradation pathways by preferentially cleaving proteins that have alanine residues. In addition, Elastase-3B is thought to play a role in the metabolism and transport of cholesterol molecules in the intestine.

REFERENCES

1. Zaks, A.S. and Bykova, A.A. 1979. Participation of the autoimmune mechanism of regulating chemical homeostasis in reactions to exposure to stress. *Patol. Fiziol. Eksp. Ter.* 4: 23-28.
2. Shirasu, Y., Takemura, K., Yoshida, H., Sato, Y., Iijima, H., Shimada, Y., Mikayama, T., Ozawa, T., Ikeda, N. and Ishida, A. 1988. Molecular cloning of complementary DNA encoding one of the human pancreatic protease E isozymes. *J. Biochem.* 104: 259-264.
3. Tani, T., Ohsumi, J., Mita, K. and Takiguchi, Y. 1988. Identification of a novel class of elastase isozyme, human pancreatic elastase III, by cDNA and genomic gene cloning. *J. Biol. Chem.* 263: 1231-1239.
4. Shimada, S., Yamaguchi, K., Takahashi, M. and Ogawa, M. 2002. Pancreatic elastase IIIA and its variants are expressed in pancreatic carcinoma cells. *Int. J. Mol. Med.* 10: 599-603.
5. Nemoda, Z. and Sahin-Tóth, M. 2006. Chymotrypsin C (caldecrin) stimulates autoactivation of human cationic trypsinogen. *J. Biol. Chem.* 281: 11879-11886.
6. Lowe, A.W., Olsen, M., Hao, Y., Lee, S.P., Taek Lee, K., Chen, X., van de Rijn, M. and Brown, P.O. 2007. Gene expression patterns in pancreatic tumors, cells and tissues. *PLoS ONE* 2: e323.
7. Teichmann, J., Mann, S.T., Stracke, H., Lange, U., Hardt, P.D., Bretzel, R.G. and Klör, H.U. 2008. Parathormone levels and Vitamin D metabolism in female patients with various grades of fecal elastase 1 deficiency. *Eur. J. Med. Res.* 13: 563-567.
8. Gao, J., Zhu, F., Lv, S., Li, Z., Ling, Z., Gong, Y., Jie, C. and Ma, L. 2010. Identification of pancreatic juice proteins as biomarkers of pancreatic cancer. *Oncol. Rep.* 23: 1683-1692.

CHROMOSOMAL LOCATION

Genetic locus: Cela3b (mouse) mapping to 4 D3.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Elastase-3 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Elastase-3 shRNA Plasmid (m): sc-144629-SH and Elastase-3 shRNA (m) Lentiviral Particles: sc-144629-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Elastase-3 siRNA (m) is recommended for the inhibition of Elastase-3 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Elastase-3 gene expression knockdown using RT-PCR Primer: Elastase-3 (m)-PR: sc-144629-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.