

Fibrocystin L siRNA (m): sc-145175

BACKGROUND

Fibrocystin is a type I membrane protein that undergoes regulated proteolysis. Many proteolytic cleavages occur on the ectodomain whereas at least one cleavage occurs on the cytoplasmic portion. Fibrocystin may participate in the mediation of intracellular calcium in the distal nephron in a manner similar to PKD1 and PKD2. Mutations in the PKHD1 gene, which encodes Fibrocystin, result in autosomal recessive polycystic kidney disease (ARPKD), a severe form of polycystic kidney disease characterized by enlarged kidneys and congenital hepatic fibrosis. A related protein, Fibrocystin L, also designated polycystic kidney and hepatic disease 1-like protein 1 or PKHD1L1, shares 41% similarity with Fibrocystin in the extracellular domain, but is not associated with ARPKD. Fibrocystin L is a large receptor protein with a signal peptide, a single transmembrane domain and a short cytoplasmic tail. It is ubiquitously expressed at low levels, with higher expression in spleen and thymus as well as in activated T cells and B lymphoblasts, suggesting a role for Fibrocystin L in the immune response.

REFERENCES

1. Ward, C.J., Hogan, M.C., Rossetti, S., Walker, D., Sneddon, T., Wang, X., Kubly, V., Cunningham, J.M., Bacallao, R., Ishibashi, M., Milliner, D.S., Torres, V.E. and Harris, P.C. 2002. The gene mutated in autosomal recessive polycystic kidney disease encodes a large, receptor-like protein. *Nat. Genet.* 30: 259-269.
2. Hogan, M.C., Griffin, M.D., Rossetti, S., Torres, V.E., Ward, C.J. and Harris, P.C. 2003. PKHD1, a homolog of the autosomal recessive polycystic kidney disease gene, encodes a receptor with inducible T lymphocyte expression. *Hum. Mol. Genet.* 12: 685-698.
3. Online Mendelian Inheritance in Man, OMIM[™]. 2004. Johns Hopkins University, Baltimore, MD. MIM Number: 607843. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Al-Bhalal, L. and Akhtar, M. 2008. Molecular basis of autosomal recessive polycystic kidney disease (ARPKD). *Adv. Anat. Pathol.* 15: 54-58.
5. Harris, P.C. and Torres, V.E. 2008. Polycystic kidney disease. *Annu. Rev. Med.* 60: 321-337.
6. Kim, I., Fu, Y., Hui, K., Moeckel, G., Mai, W., Li, C., Liang, D., Zhao, P., Ma, J., Chen, X.Z., George, A.L., Coffey, R.J., Feng, Z.P. and Wu, G. 2008. Fibrocystin/polyductin modulates renal tubular formation by regulating polycystin-2 expression and function. *J. Am. Soc. Nephrol.* 19: 455-468.
7. Kim, I., Li, C., Liang, D., Chen, X.Z., Coffey, R.J., Ma, J., Zhao, P. and Wu, G. 2008. Polycystin-2 expression is regulated by a PC2-binding domain in the intracellular portion of fibrocystin. *J. Biol. Chem.* 283: 31559-31566.
8. Fischer, D.C., Jacoby, U., Pape, L., Ward, C.J., Kuwertz-Broeking, E., Renken, C., Nizze, H., Querfeld, U., Rudolph, B., Mueller-Wiefel, D.E., Bergmann, C. and Haffner, D. 2009. Activation of the AKT/mTOR pathway in autosomal recessive polycystic kidney disease (ARPKD). *Nephrol. Dial. Transplant.* 24: 1819-1827.

CHROMOSOMAL LOCATION

Genetic locus: Pkhd1l1 (mouse) mapping to 15 B3.2.

PRODUCT

Fibrocystin L siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Fibrocystin L shRNA Plasmid (m): sc-145175-SH and Fibrocystin L shRNA (m) Lentiviral Particles: sc-145175-V as alternate gene silencing products.

For independent verification of Fibrocystin L (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-145175A, sc-145175B and sc-145175C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Fibrocystin L siRNA (m) is recommended for the inhibition of Fibrocystin L expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Fibrocystin L gene expression knockdown using RT-PCR Primer: Fibrocystin L (m)-PR: sc-145175-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.