

GalNAc-T4 siRNA (m): sc-145312

BACKGROUND

The UDP-N-acetyl- α -D-galactosamine:polypeptide N-acetylgalactosaminyltransferase (GalNAc-T) family of enzymes are substrate-specific proteins that catalyze the transfer of GalNAc (N-acetylgalactosamine) to serine and threonine residues onto various proteins, thereby initiating mucin-type O-linked glycosylation in the Golgi apparatus. GalNAc-T4 (polypeptide N-acetylgalactosaminyltransferase 4), also known as UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 4, is a 578 amino acid protein that displays high enzymatic activity toward Muc7, Muc2 and EA2. GalNAc-T4 also glycosylates PSGL-1 on threonine 57. Its N-terminal domain is involved in substrate binding and manganese coordination, while the C-terminal domain is involved in UDP-Gal binding and catalytic reaction. Though ubiquitously expressed, GalNAc-T4 is found at highest levels in mucosal cells.

REFERENCES

1. Bennett, E.P., Hassan, H. and Clausen, H. 1996. cDNA cloning and expression of a novel human UDP-N-acetyl- α -D-galactosamine. Polypeptide N-acetylgalactosaminyltransferase, GalNAc-T3. *J. Biol. Chem.* 271: 17006-17012.
2. Bennett, E.P., Hassan, H., Mandel, U., Mirgorodskaya, E., Roepstorff, P., Burchell, J., Taylor-Papadimitriou, J., Hollingsworth, M.A., Merckx, G., van Kessel, A.G., Eiberg, H., Steffensen, R. and Clausen, H. 1998. Cloning of a human UDP-N-acetyl- α -D-Galactosamine:polypeptide N-acetylgalactosaminyltransferase that complements other GalNAc-transferases in complete O-glycosylation of the MUC1 tandem repeat. *J. Biol. Chem.* 273: 30472-30481.
3. Online Mendelian Inheritance in Man, OMIM™. 1999. Johns Hopkins University, Baltimore, MD. MIM Number: 603565. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Hassan, H., Reis, C.A., Bennett, E.P., Mirgorodskaya, E., Roepstorff, P., Hollingsworth, M.A., Burchell, J., Taylor-Papadimitriou, J. and Clausen, H. 2000. The lectin domain of UDP-N-acetyl-D-galactosamine: polypeptide N-acetylgalactosaminyltransferase-T4 directs its glycopeptide specificities. *J. Biol. Chem.* 275: 38197-38205.
5. Sommer, P., Napier, H.R., Hogan, B.L. and Kidson, S.H. 2006. Identification of TGF β 1i4 as a downstream target of Foxc1. *Dev. Growth Differ.* 48: 297-308.
6. Wright, P.K., May, F.E., Darby, S., Saif, R., Lennard, T.W. and Westley, B.R. 2009. Estrogen regulates vesicle trafficking gene expression in EFF-3, EFM-19 and MCF-7 breast cancer cells. *Int. J. Clin. Exp. Pathol.* 2: 463-475.
7. O'Halloran, A.M., Patterson, C.C., Horan, P., Maree, A., Curtin, R., Stanton, A., McKeown, P.P. and Shields, D.C. 2009. Genetic polymorphisms in platelet-related proteins and coronary artery disease: investigation of candidate genes, including N-acetylgalactosaminyltransferase 4 (GALNT4) and sulphotransferase 1A1/2 (SULT1A1/2). *J. Thromb. Thrombolysis* 27: 175-184.
8. Lu, L., Niu, B., Zhao, J., Liu, L., Lu, W.C., Liu, X.J., Li, Y.X. and Cai, Y.D. 2009. GalNAc-transferase specificity prediction based on feature selection method. *Peptides* 30: 359-364.

CHROMOSOMAL LOCATION

Genetic locus: Galnt4 (mouse) mapping to 10 C3.

PRODUCT

GalNAc-T4 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GalNAc-T4 shRNA Plasmid (m): sc-145312-SH and GalNAc-T4 shRNA (m) Lentiviral Particles: sc-145312-V as alternate gene silencing products.

For independent verification of GalNAc-T4 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-145312A, sc-145312B and sc-145312C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GalNAc-T4 siRNA (m) is recommended for the inhibition of GalNAc-T4 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GalNAc-T4 gene expression knockdown using RT-PCR Primer: GalNAc-T4 (m)-PR: sc-145312-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.