

GIGYF2 siRNA (m): sc-145397

BACKGROUND

GIGYF2 (GRB10 interacting GYF protein 2), also known as GYF2, PERQ2, PERQ3, PARK11 or TNRC15, is a 1,299 amino acid protein that may be involved in the regulation of tyrosine kinase receptor signaling, including IGF-I and Insulin receptors. Belonging to the PERQ family of proteins, GIGYF2 contains long stretches of glutamine and glutamic acid residues. Mutations in the gene encoding GIGYF2 are the cause of Parkinson disease type 11 (PARK11), which is characterized by bradykinesia, resting tremor, muscular rigidity and postural instability. Parkinson's disease involves the loss of dopaminergic neurons in the substantia nigra and the presence of Lewy bodies (intraneuronal accumulations of aggregated proteins), in surviving neurons in various areas of the brain. PARK11 may show age-dependent or reduced penetrance. GIGYF2 exists as two alternatively spliced isoforms.

REFERENCES

1. Lautier, C., Goldwurm, S., Dürr, A., Giovannone, B., Tsiaras, W.G., Pezzoli, G., Brice, A. and Smith, R.J. 2008. Mutations in the GIGYF2 (TNRC15) gene at the PARK11 locus in familial Parkinson disease. *Am. J. Hum. Genet.* 82: 822-833.
2. Bonifati, V. 2009. Is GIGYF2 the defective gene at the PARK11 locus? *Curr. Neurol. Neurosci. Rep.* 9: 185-187.
3. Bras, J., Simón-Sánchez, J., Federoff, M., Morgado, A., Januario, C., Ribeiro, M., Cunha, L., Oliveira, C. and Singleton, A.B. 2009. Lack of replication of association between GIGYF2 variants and Parkinson disease. *Hum. Mol. Genet.* 18: 341-346.
4. Sutherland, G.T., Siebert, G.A., Newman, J.R., Silburn, P.A., Boyle, R.S., O'Sullivan, J.D. and Mellick, G.D. 2009. Haplotype analysis of the PARK 11 gene, GIGYF2, in sporadic Parkinson's disease. *Mov. Disord.* 24: 449-452.
5. Lesage, S., Condroyer, C., Lohman, E., Troiano, A., Tison, F., Viallet, F., Damier, P., Tranchant, C., Vidhaillet, M., Ouvrard-Hernandez, A.M., Dürr, A. and Brice, A. 2009. Follow-up study of the GIGYF2 gene in French families with Parkinson's disease. *Neurobiol. Aging* 31: 1069-1071.

CHROMOSOMAL LOCATION

Genetic locus: Gygf2 (mouse) mapping to 1 D.

PRODUCT

GIGYF2 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GIGYF2 shRNA Plasmid (m): sc-145397-SH and GIGYF2 shRNA (m) Lentiviral Particles: sc-145397-V as alternate gene silencing products.

For independent verification of GIGYF2 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-145397A, sc-145397B and sc-145397C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GIGYF2 siRNA (m) is recommended for the inhibition of GIGYF2 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

GIGYF2 (A-12): sc-393918 is recommended as a control antibody for monitoring of GIGYF2 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GIGYF2 gene expression knockdown using RT-PCR Primer: GIGYF2 (m)-PR: sc-145397-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.