



Glut7 siRNA (m): sc-145450

BACKGROUND

Glucose serves as the major energy substrate of mammalian cells and is fundamental to metabolism. Glucose passage across cell membranes is mediated by a family of transporters termed glucose transporters, or Gluts, which are characterized by the presence of 12 membrane-spanning helices. The Glut family is divided into three subfamilies: class I (previously known as glucose transporters), which includes Glut1, Glut2, Glut3 and Glut4; class II (previously known as fructose transporters), which includes Glut5, Glut7, Glut9 and Glut11; and class III, which includes Glut6, Glut8, Glut10, Glut12 and the myo-inositol transporter HMIT1. Glut7, also known as SLC2A7 (solute carrier family 2 (facilitated glucose transporter), member 7), is a 513 amino acid multi-pass membrane protein belonging to the sugar transporter family, glucose transporter subfamily and major facilitator superfamily.

REFERENCES

1. Bell, G.I., Kayano, T., Buse, J.B., Burant, C.F., Takeda, J., Lin, D., Fukumoto, H. and Seino, S. 1990. Molecular biology of mammalian glucose transporters. *Diabetes Care* 13: 198-208.
2. Mueckler, M. 1994. Facilitative glucose transporters. *Eur. J. Biochem.* 219: 713-725.
3. Joost, H.G. and Thorens, B. 2001. The extended Glut-family of sugar/polyol transport facilitators: nomenclature, sequence characteristics, and potential function of its novel members (review). *Mol. Membr. Biol.* 18: 247-256.
4. Li, Q., Manolescu, A., Ritzel, M., Yao, S., Slugoski, M., Young, J.D., Chen, X.Z. and Cheeseman, C.I. 2004. Cloning and functional characterization of the human Glut7 isoform SLC2A7 from the small intestine. *Am. J. Physiol. Gastrointest. Liver Physiol.* 287: G236-G242.
5. Manolescu, A., Salas-Burgos, A.M., Fischbarg, J. and Cheeseman, C.I. 2005. Identification of a hydrophobic residue as a key determinant of fructose transport by the facilitative hexose transporter SLC2A7 (Glut7). *J. Biol. Chem.* 280: 42978-42983.
6. Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 610371. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. Manolescu, A.R., Witkowska, K., Kinnaird, A., Cessford, T. and Cheeseman, C. 2007. Facilitated hexose transporters: new perspectives on form and function. *Physiology* 22: 234-240.
8. Cheeseman, C. 2008. Glut7: a new intestinal facilitated hexose transporter. *Am. J. Physiol. Endocrinol. Metab.* 295: E238-E241.

CHROMOSOMAL LOCATION

Genetic locus: Slc2a7 (mouse) mapping to 4 E2.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Glut7 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Glut7 shRNA Plasmid (m): sc-145450-SH and Glut7 shRNA (m) Lentiviral Particles: sc-145450-V as alternate gene silencing products.

For independent verification of Glut7 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-145450A, sc-145450B and sc-145450C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Glut7 siRNA (m) is recommended for the inhibition of Glut7 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Glut7 gene expression knockdown using RT-PCR Primer: Glut7 (m)-PR: sc-145450-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.