

# FAM170A siRNA (m): sc-145634

## BACKGROUND

FAM170A (family with sequence similarity 170, member A), also known as ZNFD (zinc finger domain-containing protein), is a 330 amino acid protein containing one C<sub>2</sub>H<sub>2</sub>-type zinc finger. Existing as four alternatively spliced isoforms, FAM170A is encoded by a gene located on human chromosome 5q23.1. Chromosome 5 consists of 181 million base pairs and comprises nearly 6% of the human genome. Chromosome 5 is associated with Cockayne syndrome through the ERCC8 gene and familial adenomatous polyposis through the adenomatous polyposis coli (APC) tumor suppressor gene. Treacher Collins syndrome is also chromosome 5-associated and is caused by insertions or deletions within the TCOF1 gene. Deletion of the p arm of chromosome 5 leads to Cri du chat syndrome, while deletion of the q arm of chromosome 5 altogether is common in therapy-related acute myelogenous leukemias and myelodysplastic syndrome.

## REFERENCES

1. Rauch, A. and Dörr, H.G. 2007. Chromosome 5q subtelomeric deletion syndrome. *Am. J. Med. Genet. C Semin. Med. Genet.* 145C: 372-376.
2. Villa, N., Redaelli, S., Borroni, C., Colombo, C., Roncaglia, N., Sala, E., Crosti, F., Cappellini, A. and Dalprà, L. 2007. Fetal trisomy 5 mosaicism: case report and literature review. *Am. J. Med. Genet. A* 143A: 2343-2346.
3. Shadduck, R.K., Latsko, J.M., Rossetti, J.M., Haq, B. and Abdulhaq, H. 2007. Recent advances in myelodysplastic syndromes. *Exp. Hematol.* 35: 137-143.
4. Falini, B., Nicoletti, I., Bolli, N., Martelli, M.P., Liso, A., Gorello, P., Mandelli, F., Mecucci, C. and Martelli, M.F. 2007. Translocations and mutations involving the nucleophosmin (NPM1) gene in lymphomas and leukemias. *Haematologica* 92: 519-532.
5. Kristoffersen, K.E. 2008. Speech and language development in cri du chat syndrome: a critical review. *Clin. Linguist. Phon.* 22: 443-457.
6. Valent, P. 2008. Revealing the pathogenesis of the 5q- syndrome. *Eur. J. Clin. Invest.* 38: 539-540.
7. Buysse, K., Crepel, A., Menten, B., Pattyn, F., Antonacci, F., Veltman, J.A., Larsen, L.A., Tümer, Z., de Klein, A., van de Laar, I., Devriendt, K., Mortier, G. and Speleman, F. 2008. Mapping of 5q35 chromosomal rearrangements within a genomically unstable region. *J. Med. Genet.* 45: 672-678.
8. Azman, B.Z., Akhir, S.M., Zilfalil, B.A. and Ankathil, R. 2008. Two cases of deletion 5p syndrome: one with paternal involvement and another with atypical presentation. *Singapore Med. J.* 49: e98-e100.

## CHROMOSOMAL LOCATION

Genetic locus: Fam170a (mouse) mapping to 18 D1.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

FAM170A siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 µM solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see FAM170A shRNA Plasmid (m): sc-145634-SH and FAM170A shRNA (m) Lentiviral Particles: sc-145634-V as alternate gene silencing products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 µl of the RNase-free water provided. Resuspension of the siRNA duplex in 330 µl of RNase-free water makes a 10 µM solution in a 10 µM Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

FAM170A siRNA (m) is recommended for the inhibition of FAM170A expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 µM in 66 µl. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor FAM170A gene expression knockdown using RT-PCR Primer: FAM170A (m)-PR: sc-145634-PR (20 µl). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.