



HATH-6 siRNA (m): sc-145899

BACKGROUND

HATH-6, also known as atonal homolog 8 (ATOH8), is a 321 amino acid putative transcription factor. Localized to the nucleus, HATH-6 is expressed in kidney, lung, liver, pancreas, heart and the endothelium of umbilical vessels. HATH-6 is thought to participate in podocyte development in the kidney and may be involved in specification and differentiation of brain neuronal cell lineages. As a transcription factor, HATH-6 contains one basic helix-loop-helix (bHLH) domain. In order to efficiently bind to DNA, HATH-6 must be bound to another bHLH domain containing protein. HATH-6 exists as two isoforms produced by alternative splicing.

REFERENCES

1. Murre, C., McCaw, P.S. and Baltimore, D. 1989. A new DNA binding and dimerization motif in immunoglobulin enhancer binding, daughterless, MyoD, and myc proteins. *Cell* 56: 777-783.
2. Kato, G.J. and Dang, C.V. 1992. Function of the c-Myc oncoprotein. *FASEB J.* 6: 3065-3072.
3. Inoue, C., Bae, S.K., Takatsuka, K., Inoue, T., Bessho, Y. and Kageyama, R. 2001. Math6, a bHLH gene expressed in the developing nervous system, regulates neuronal versus glial differentiation. *Genes Cells* 6: 977-986.
4. Ledent, V. and Vervoort, M. 2001. The basic helix-loop-helix protein family: comparative genomics and phylogenetic analysis. *Genome Res.* 11: 754-770.
5. Ledent, V., Paquet, O. and Vervoort, M. 2002. Phylogenetic analysis of the human basic helix-loop-helix proteins. *Genome Biol.* 3: RESEARCH0030.
6. Wasserman, S.M., Mehraban, F., Komuves, L.G., Yang, R.B., Tomlinson, J.E., Zhang, Y., Spriggs, F. and Topper, J.N. 2002. Gene expression profile of human endothelial cells exposed to sustained fluid shear stress. *Physiol. Genomics* 12: 13-23.
7. Ross, M.D., Martinka, S., Mukherjee, A., Sedor, J.R., Vinson, C. and Bruggeman, L.A. 2006. Math6 expression during kidney development and altered expression in a mouse model of glomerulosclerosis. *Dev. Dyn.* 235: 3102-3109.
8. Kautz, L., Meynard, D., Monnier, A., Darnaud, V., Bouvet, R., Wang, R.H., Deng, C., Vaulont, S., Mosser, J., Coppin, H. and Roth, M.P. 2008. Iron regulates phosphorylation of Smad1/5/8 and gene expression of Bmp6, Smad7, Id1, and Atoh8 in the mouse liver. *Blood* 112: 1503-1509.

CHROMOSOMAL LOCATION

Genetic locus: Atoh8 (mouse) mapping to 6 C1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

HATH-6 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see HATH-6 shRNA Plasmid (m): sc-145899-SH and HATH-6 shRNA (m) Lentiviral Particles: sc-145899-V as alternate gene silencing products.

For independent verification of HATH-6 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-145899A and sc-145899B.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

HATH-6 siRNA (m) is recommended for the inhibition of HATH-6 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor HATH-6 gene expression knockdown using RT-PCR Primer: HATH-6 (m)-PR: sc-145899-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.