

Histone cluster 1 H3E siRNA (m): sc-146007

BACKGROUND

Eukaryotic histones are basic and water soluble nuclear proteins that form hetero-octameric nucleosome particles by wrapping 146 base pairs of DNA in a left-handed super-helical turn sequentially to form chromosomal fiber. Two molecules of each of the four core Histones (H2A, H2B, H3, and H4) form the octamer, which is comprised of two H2A-H2B dimers and two H3-H4 dimers, forming two nearly symmetrical halves by tertiary structure. Histones are subject to posttranslational modification by enzymes primarily on their N-terminal tails, but also in their globular domains. Histone cluster 1 H3E is a member of the H3 Histone family and is encoded by a gene located in the large histone gene cluster on chromosome 6p22.2. Histone cluster 1 H3E undergoes posttranslational modifications, including the phosphorylation of Ser 10.

REFERENCES

1. Albig, W., Kardalidou, E., Drabent, B., Zimmer, A. and Doenecke, D. 1991. Isolation and characterization of two human H1 Histone genes within clusters of core histone genes. *Genomics* 10: 940-948.
2. Albig, W., Kioschis, P., Poustka, A., Meergans, K. and Doenecke, D. 1997. Human histone gene organization: nonregular arrangement within a large cluster. *Genomics* 40: 314-322.
3. Albig, W. and Doenecke, D. 1997. The human histone gene cluster at the D6S105 locus. *Hum. Genet.* 101: 284-294.
4. Goto, H., Tomono, Y., Ajiro, K., Kosako, H., Fujita, M., Sakurai, M., Okawa, K., Iwamatsu, A., Okigaki, T., Takahashi, T. and Inagaki, M. 1999. Identification of a novel phosphorylation site on Histone H3 coupled with mitotic chromosome condensation. *J. Biol. Chem.* 274: 25543-25549.
5. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 602813. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Wang, H., Zhai, L., Xu, J., Joo, H.Y., Jackson, S., Erdjument-Bromage, H., Tempst, P., Xiong, Y. and Zhang, Y. 2006. Histone H3 and H4 ubiquitylation by the CUL4-DDB-ROC1 ubiquitin ligase facilitates cellular response to DNA damage. *Mol. Cell* 22: 383-394.
7. Beck, H.C., Nielsen, E.C., Matthiesen, R., Jensen, L.H., Sehested, M., Finn, P., Grauslund, M., Hansen, A.M. and Jensen, O.N. 2006. Quantitative proteomic analysis of post-translational modifications of human histones. *Mol. Cell. Proteomics* 5: 1314-1325.
8. Morris, S.A., Rao, B., Garcia, B.A., Hake, S.B., Diaz, R.L., Shabanowitz, J., Hunt, D.F., Allis, C.D., Lieb, J.D. and Strahl, B.D. 2007. Identification of Histone H3 lysine 36 acetylation as a highly conserved histone modification. *J. Biol. Chem.* 282: 7632-7640.
9. Garcia, B.A., Hake, S.B., Diaz, R.L., Kauer, M., Morris, S.A., Recht, J., Shabanowitz, J., Mishra, N., Strahl, B.D., Allis, C.D. and Hunt, D.F. 2007. Organismal differences in post-translational modifications in Histones H3 and H4. *J. Biol. Chem.* 282: 7641-7655.

CHROMOSOMAL LOCATION

Genetic locus: Hist1h3e (mouse) mapping to 13 A3.1.

PRODUCT

Histone cluster 1 H3E siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Histone cluster 1 H3E shRNA Plasmid (m): sc-146007-SH and Histone cluster 1 H3E shRNA (m) Lentiviral Particles: sc-146007-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Histone cluster 1 H3E siRNA (m) is recommended for the inhibition of Histone cluster 1 H3E expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.