

IL-1F8 siRNA (m): sc-146211

BACKGROUND

IL-1 (Interleukin-1) is a cytokine responsible for initiating a variety of activities through the activation of transcription factors, NF κ B and AP-1, thereby promoting host response to injury or infection. The IL-1 superfamily is comprised of several ligands and receptors. IL-1F8, also known as IL-1h, or interleukin-1 homolog 2 (IL-1H2), is a secreted ligand belonging to this superfamily. IL-1F8 is highly expressed in epithelial cells but is also found in skeletal muscle and glial cells. IL-1F8 activates the IL-1Rrp2 and IL-1RAcP-dependent pathways leading to MAPKs, NF κ B activation and stimulation of IL-6 and IL-8 production. In addition, IL-1F8 may play a role in the pathogenesis of rheumatoid arthritis (RA). Two isoforms exist for this protein, isoform 1 and isoform 2. These isoforms differ from one another in their amino acid sequences between residues 88 and 164.

REFERENCES

1. Kumar, S., McDonnell, P.C., Lehr, R., Tierney, L., Tzimas, M.N., Griswold, D.E., Capper, E.A., Tal-Singer, R., Wells, G.I., Doyle, M.L. and Young, P.R. 2000. Identification and initial characterization of four novel members of the interleukin-1 family. *J. Biol. Chem.* 275: 10308-10314.
2. Gao, W., Kumar, S., Lotze, M.T., Hanning, C., Robbins, P.D. and Gambotto, A. 2002. Innate immunity mediated by the cytokine IL-1 homologue 4 (IL-1H4/IL-1F7) induces IL-12-dependent adaptive and profound antitumor immunity. *J. Immunol.* 170: 107-113.
3. Jensen, L.E. and Whitehead, A.S. 2004. The 3' untranslated region of the membrane-bound IL-1 accessory protein mRNA confers tissue-specific destabilization. *J. Immunol.* 173: 6248-6258.
4. Timms, A.E., Crane, A.M., Sims, A.M., Cordell, H.J., Bradbury, L.A., Abbott, A., Coyne, M.R., Beynon, O., Herzberg, I., Duff, G.W., Calin, A., Cardon, L.R., Wordsworth, B.P. and Brown, M.A. 2004. The interleukin 1 gene cluster contains a major susceptibility locus for ankylosing spondylitis. *Am. J. Hum. Genet.* 75: 587-595.
5. Towne, J.E., Garka, K.E., Renshaw, B.R., Virca, G.D. and Sims, J.E. 2004. Interleukin (IL)-1F6, IL-1F8, and IL-1F9 signal through IL-1Rrp2 and IL-1RAcP to activate the pathway leading to NF κ B and MAPKs. *J. Biol. Chem.* 279: 13677-13688.
6. Wang, P., Meinhardt, B., Andre, R., Renshaw, B.R., Kimber, I., Rothwell, N.J. and Pinteaux, E. 2005. The interleukin-1-related cytokine IL-1F8 is expressed in glial cells, but fails to induce IL-1 β signalling responses. *Cytokine* 29: 245-250.
7. Magne, D., Palmer, G., Barton, J.L., Mezin, F., Talabot-Ayer, D., Bas, S., Duffy, T., Noger, M., Guerne, P.A., Nicklin, M.J. and Gabay, C. 2006. The new IL-1 family member IL-1F8 stimulates production of inflammatory mediators by synovial fibroblasts and articular chondrocytes. *Arthritis Res. Ther.* 8: R80.
8. Barksby, H.E., Lea, S.R., Preshaw, P.M. and Taylor, J.J. 2007. The expanding family of interleukin-1 cytokines and their role in destructive inflammatory disorders. *Clin. Exp. Immunol.* 149: 217-225.
9. Knedla, A., Neumann, E. and Müller-Ladner, U. 2007. Developments in the synovial biology field 2006. *Arthritis Res. Ther.* 9: 209.

CHROMOSOMAL LOCATION

Genetic locus: IL1f8 (mouse) mapping to 2 A3.

PRODUCT

IL-1F8 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see IL-1F8 shRNA Plasmid (m): sc-146211-SH and IL-1F8 shRNA (m) Lentiviral Particles: sc-146211-V as alternate gene silencing products.

For independent verification of IL-1F8 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-146211A and sc-146211B.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

IL-1F8 siRNA (m) is recommended for the inhibition of IL-1F8 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor IL-1F8 gene expression knockdown using RT-PCR Primer: IL-1F8 (m)-PR: sc-146211-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.