

## IL-34 siRNA (m): sc-146220

### BACKGROUND

The interleukins (ILs) are a broad family of well characterized cytokines, primarily of hematopoietic cell origin. They are secreted by immune cells (mainly macrophages, B cells or T cells) that regulate a wide range of immune system functions. The specific functions of different interleukins vary from the regulation of inflammatory and immune responses to the regulation of other interleukins. IL-34 (interleukin-34) is a 235 amino acid secreted mouse protein that belongs to the interleukin family of cytokines. Existing as a homodimer, IL-34 functions to promote the viability and differentiation of macrophages and monocytes and may also act as a ligand for c-Fms/CSF-1R, a colony-stimulating factor. Multiple alternatively spliced isoforms of IL-34 exist, all of which are encoded by a gene that maps to mouse chromosome 8.

### REFERENCES

1. Smith, K.A., Lachman, L.B., Oppenheim, J.J. and Favata, M.F. 1980. The functional relationship of the interleukins. *J. Exp. Med.* 151: 1551-1556.
2. Cockayne, D.A., Abrams, J.S. and Nienhuis, A.W. 1991. Antisense RNA inhibition of hematopoietic growth factor production. *Growth Factors* 5: 171-181.
3. Sander, B., Höiden, I., Andersson, U., Möller, E. and Abrams, J.S. 1993. Similar frequencies and kinetics of cytokine producing cells in murine peripheral blood and spleen. Cytokine detection by immunoassay and intracellular immunostaining. *J. Immunol. Methods* 166: 201-214.
4. Moldenhauer, A., Genter, G., Lun, A., Bal, G., Kiesewetter, H. and Salama, A. 2008. Hematopoietic progenitor cells and interleukin-stimulated endothelium: expansion and differentiation of myeloid precursors. *BMC Immunol.* 9: 56.
5. Lin, H., Lee, E., Hestir, K., Leo, C., Huang, M., Bosch, E., Halenbeck, R., Wu, G., Zhou, A., Behrens, D., Hollenbaugh, D., Linnemann, T., Qin, M., Wong, J., Chu, K., Doberstein, S.K. and Williams, L.T. 2008. Discovery of a cytokine and its receptor by functional screening of the extracellular proteome. *Science* 320: 807-811.
6. Liptrott, N.J., Penny, M., Bray, P.G., Sathish, J., Khoo, S.H., Back, D.J. and Owen, A. 2009. The impact of cytokines on the expression of drug transporters, cytochrome P450 enzymes and chemokine receptors in human PBMC. *Br. J. Pharmacol.* 156: 497-508.
7. Cheng, M., Charoudeh, H.N., Brodin, P., Tang, Y., Lakshmikanth, T., Höglund, P., Jacobsen, S.E. and Sitnicka, E. 2009. Distinct and overlapping patterns of cytokine regulation of thymic and bone marrow-derived NK cell development. *J. Immunol.* 182: 1460-1468.

### CHROMOSOMAL LOCATION

Genetic locus: Il34 (mouse) mapping to 8 E1.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

### PRODUCT

IL-34 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see IL-34 shRNA Plasmid (m): sc-146220-SH and IL-34 shRNA (m) Lentiviral Particles: sc-146220-V as alternate gene silencing products.

For independent verification of IL-34 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-146220A and sc-146220B.

### STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

### APPLICATIONS

IL-34 siRNA (m) is recommended for the inhibition of IL-34 expression in mouse cells.

### SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

### RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor IL-34 gene expression knockdown using RT-PCR Primer: IL-34 (m)-PR: sc-146220-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.