

ISG20L1 siRNA (m): sc-146296

BACKGROUND

ISG20L1 (apoptosis enhancing nuclease), also known as AEN or pp12744, is a 325 amino acid protein that contains one exonuclease domain and shares significant homology with various exonucleases. Encoded by a gene that maps to human chromosome 15q26.1, ISG20L1 localizes predominantly to nucleolus where it translocates to nucleoplasm upon apoptosis induction. Up-regulated by p53 in response to ionizing radiation and DNA-damaging agents, ISG20L1 requires phosphorylation of p53 at Ser 15 for effective induction. ISG20L1 digests double-stranded DNA to form single-stranded DNA and amplifies DNA damage signals, thereby mediating and enhancing p53-induced apoptosis. Existing as two alternatively spliced isoforms, ISG20L1 exhibits exonuclease activity against single- and double-stranded DNA and RNA and also participates in nucleic acid binding.

REFERENCES

1. Lee, J.H., Koh, Y.A., Cho, C.K., Lee, S.J., Lee, Y.S. and Bae, S. 2005. Identification of a novel ionizing radiation-induced nuclease, AEN, and its functional characterization in apoptosis. *Biochem. Biophys. Res. Commun.* 337: 39-47.
2. Turtoi, A., Brown, I., Oskamp, D. and Schneeweiss, F.H. 2008. Early gene expression in human lymphocytes after γ -irradiation—a genetic pattern with potential for biodosimetry. *Int. J. Radiat. Biol.* 84: 375-387.
3. Coute, Y., Kindbeiter, K., Belin, S., Dieckmann, R., Duret, L., Bezin, L., Sanchez, J.C. and Diaz, J.J. 2008. ISG20L2, a novel vertebrate nucleolar exoribonuclease involved in ribosome biogenesis. *Mol. Cell. Proteomics* 7: 546-559.
4. Kawase, T., Ichikawa, H., Ohta, T., Nozaki, N., Tashiro, F., Ohki, R. and Taya, Y. 2008. p53 target gene AEN is a nuclear exonuclease required for p53-dependent apoptosis. *Oncogene* 27: 3797-3810.
5. Zauli, G., di Iasio, M.G., Secchiero, P., Dal Bo, M., Marconi, D., Bomben, R., Del Poeta, G. and Gattei, V. 2009. Exposure of B cell chronic lymphocytic leukemia (B-CLL) cells to nutlin-3 induces a characteristic gene expression profile, which correlates with nutlin-3-mediated cytotoxicity. *Curr. Cancer Drug Targets* 9: 510-518.
6. Turtoi, A. and Schneeweiss, F.H. 2009. Effect of (211)At α -particle irradiation on expression of selected radiation responsive genes in human lymphocytes. *Int. J. Radiat. Biol.* 85: 403-412.
7. Marsh, B.J., Stevens, S.L., Hunter, B. and Stenzel-Poore, M.P. 2009. Inflammation and the emerging role of the Toll-like receptor system in acute brain ischemia. *Stroke* 40: S34-S37.
8. Haralambieva, I.H., Dhiman, N., Ovsyannikova, I.G., Vierkant, R.A., Pankratz, V.S., Jacobson, R.M. and Poland, G.A. 2010. 2'-5'-oligoadenylate synthetase single-nucleotide polymorphisms and haplotypes are associated with variations in immune responses to rubella vaccine. *Hum. Immunol.* 71: 383-391.
9. Eby, K.G., Rosenbluth, J.M., Mays, D.J., Marshall, C.B., Barton, C.E., Sinha, S., Johnson, K.N., Tang, L. and Pietenpol, J.A. 2010. ISG20L1 is a p53 family target gene that modulates genotoxic stress-induced autophagy. *Mol. Cancer* 9: 95.

CHROMOSOMAL LOCATION

Genetic locus: Aen (mouse) mapping to 7 D3.

PRODUCT

ISG20L1 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ISG20L1 shRNA Plasmid (m): sc-146296-SH and ISG20L1 shRNA (m) Lentiviral Particles: sc-146296-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ISG20L1 siRNA (m) is recommended for the inhibition of ISG20L1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ISG20L1 gene expression knockdown using RT-PCR Primer: ISG20L1 (m)-PR: sc-146296-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.