

# KIF18A siRNA (m): sc-146467

## BACKGROUND

The kinesins constitute a large family of microtubule-dependent motor proteins, which are responsible for the distribution of numerous organelles, vesicles and macromolecular complexes throughout the cell. Individual kinesin members play crucial roles in cell division, intracellular transport and membrane trafficking events, including endocytosis and transcytosis. KIF18A (kinesin family member 18A), also designated MS-KIF18A, is a 898 amino acid protein that contains one kinesin-motor domain, which is responsible for the ATP-dependent movement of KIF18A across microtubules. KIF18A plays a role in chromosome congression by reducing the amplitude of pre-anaphase oscillations and slowing poleward movement during anaphase, thus suppressing chromosome movements. Expression of KIF18A is induced by estrogen.

## REFERENCES

1. Luboshits, G. and Benayahu, D. 2005. MS-KIF18A, new kinesin; structure and cellular expression. *Gene* 351: 19-28.
2. Mayr, M.I., Hümmer, S., Bormann, J., Grüner, T., Adio, S., Woehlke, G. and Mayer, T.U. 2007. The human kinesin Kif18A is a motile microtubule depolymerase essential for chromosome congression. *Curr. Biol.* 17: 488-498.
3. Luboshits, G. and Benayahu, D. 2007. MS-KIF18A, a kinesin, is associated with estrogen receptor. *J. Cell. Biochem.* 100: 693-702.
4. Stumpff, J., von Dassow, G., Wagenbach, M., Asbury, C. and Wordeman, L. 2008. The kinesin-8 motor Kif18A suppresses kinetochore movements to control mitotic chromosome alignment. *Dev. Cell* 14: 252-262.
5. Zusev, M. and Benayahu, D. 2008. New insights on cellular distribution, microtubule interactions and post-translational modifications of MS-KIF18A. *J. Cell. Physiol.* 217: 618-625.
6. Gardner, M.K., Odde, D.J. and Bloom, K. 2008. Kinesin-8 molecular motors: putting the brakes on chromosome oscillations. *Trends Cell Biol.* 18: 307-310.
7. Zusev, M. and Benayahu, D. 2009. The regulation of MS-KIF18A expression and cross talk with estrogen receptor. *PLoS ONE* 4: e6407.
8. Du, Y., English, C.A. and Ohi, R. 2010. The kinesin-8 Kif18A dampens microtubule plus-end dynamics. *Curr. Biol.* 20: 374-380.

## CHROMOSOMAL LOCATION

Genetic locus: Kif18a (mouse) mapping to 2 E3.

## PRODUCT

KIF18A siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see KIF18A shRNA Plasmid (m): sc-146467-SH and KIF18A shRNA (m) Lentiviral Particles: sc-146467-V as alternate gene silencing products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

KIF18A siRNA (m) is recommended for the inhibition of KIF18A expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

KIF18A (C-12): sc-390600 is recommended as a control antibody for monitoring of KIF18A gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor KIF18A gene expression knockdown using RT-PCR Primer: KIF18A (m)-PR: sc-146467-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.