



LNK3 siRNA (m): sc-146770

BACKGROUND

The eukaryotic PDZ domain is a multifunctional protein-protein interacting motif that is found in a variety of proteins and is involved in both the clustering of signaling molecules and the organization of protein networks. LNK3 (ligand of Numb protein X 3), also known as PDZRN3 (PDZ domain containing ring finger 3) or SEMCAP3 (semaphorin cytoplasmic domain-associated protein 3), is a 1,066 amino acid protein that contains one TRAF-type zinc finger, one RING-type zinc finger and two PDZ domains. Expressed in a variety of tissues, LNK3 interacts with neuroligin 1 and ephrin-B2 and may exhibit tumor suppressive activity in ovarian serous papillary tumors. Multiple isoforms of LNK3 exist due to alternative splicing events.

REFERENCES

1. Kikuno, R., Nagase, T., Ishikawa, K., Hirose, M., Miyajima, N., Tanaka, A., Kotani, H., Nomura, N. and Ohara, O. 1999. Prediction of the coding sequences of unidentified human genes. XIV. The complete sequences of 100 new cDNA clones from brain which code for large proteins *in vitro*. DNA Res. 6: 197-205.
2. Santin, A.D., Zhan, F., Bellone, S., Palmieri, M., Cane, S., Bignotti, E., Anfossi, S., Gokden, M., Dunn, D., Roman, J.J., O'Brien, T.J., Tian, E., Cannon, M.J., Shaughnessy, J. and Pecorelli, S. 2004. Gene expression profiles in primary ovarian serous papillary tumors and normal ovarian epithelium: identification of candidate molecular markers for ovarian cancer diagnosis and therapy. Int. J. Cancer 112: 14-25.
3. Katoh, M. and Katoh, M. 2004. Identification and characterization of PDZRN3 and PDZRN4 genes in silico. Int. J. Mol. Med. 13: 607-613.
4. Meyer, G., Varoqueaux, F., Neeb, A., Oeschlies, M. and Brose, N. 2004. The complexity of PDZ domain-mediated interactions at glutamatergic synapses: a case study on neuroligin. Neuropharmacology 47: 724-733.
5. Online Mendelian Inheritance in Man, OMIM™. 2005. Johns Hopkins University, Baltimore, MD. MIM Number: 609729. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Ko, J.A., Kimura, Y., Matsuura, K., Yamamoto, H., Gondo, T. and Inui, M. 2006. PDZRN3 (LNK3, SEMCAP3) is required for the differentiation of C2C12 myoblasts into myotubes. J. Cell Sci 119: 5106-5113.

CHROMOSOMAL LOCATION

Genetic locus: Pdzrn3 (mouse) mapping to 6 D3.

PRODUCT

LNK3 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see LNK3 shRNA Plasmid (m): sc-146770-SH and LNK3 shRNA (m) Lentiviral Particles: sc-146770-V as alternate gene silencing products.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

LNK3 siRNA (m) is recommended for the inhibition of LNK3 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor LNK3 gene expression knockdown using RT-PCR Primer: LNK3 (m)-PR: sc-146770-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.