

LPAAT- η siRNA (m): sc-149019

BACKGROUND

Phosphatidic acid and lysophosphatidic acid are phospholipids involved in lipid biosynthesis and signal transduction. LPAAT- η , also known as lysophospholipid acyltransferase LPCAT4, AGPAT7 (1-acylglycerol-3-phosphate O-acyltransferase 7), AYTL3 (acyltransferase-like 3) or LPEAT2 (lysophosphatidylethanolamine acyltransferase 2), is a 524 amino acid protein belonging to the 1-acyl-sn-glycerol-3-phosphate acyltransferase family. LPAAT- η displays acyl-CoA-dependent lysophospholipid acyltransferase activity, with lysophospholipids as its substrates. For example, LPAAT- η converts lysophosphatidylethanolamine to phosphatidylethanolamine and lysophosphatidylcholine to phosphatidylcholine, respectively. In contrast, LPAAT- η has no lysophosphatidylinositol, glycerol-3-phosphate, diacylglycerol or lysophosphatidic acid acyltransferase activity. LPAAT- η also prefers long chain acyl-CoAs (C16, C18) as acyl donors. Localized to the endoplasmic reticulum membrane, LPAAT- η is widely expressed with predominant levels in brain. Two isoforms of LPAAT- η are produced by alternative splicing events.

REFERENCES

1. Siegel, H.I. and Rosenblatt, J.S. 1975. Latency and duration of estrogen induction of maternal behavior in hysterectomized-ovariectomized virgin rats: effects of pup stimulation. *Physiol. Behav.* 14: 473-476.
2. West, J., Tompkins, C.K., Balantac, N., Nudelman, E., Meengs, B., White, T., Bursten, S., Coleman, J., Kumar, A., Singer, J.W. and Leung, D.W. 1997. Cloning and expression of two human lysophosphatidic acid acyltransferase cDNAs that enhance cytokine-induced signaling responses in cells. *DNA Cell Biol.* 16: 691-701.
3. Aguado, B. and Campbell, R.D. 1998. Characterization of a human lysophosphatidic acid acyltransferase that is encoded by a gene located in the class III region of the human major histocompatibility complex. *J. Biol. Chem.* 273: 4096-4105.
4. Bursten, S.L. 1998. Interaction of lipopolysaccharide with a mammalian lyso-phosphatidate acyltransferase (LPAAT) transfected into *E. coli*, and effect of Iisofylline on LPAAT transfected into mammalian cells. *Prog. Clin. Biol. Res.* 397: 345-356.
5. Eberhardt, C., Gray, P.W. and Tjoelker, L.W. 1999. cDNA cloning, expression and chromosomal localization of two human lysophosphatidic acid acyltransferases. *Adv. Exp. Med. Biol.* 469: 351-356.
6. Ye, G.M., Chen, C., Huang, S., Han, D.D., Guo, J.H., Wan, B. and Yu, L. 2005. Cloning and characterization a novel human 1-acyl-sn-glycerol-3-phosphate acyltransferase gene AGPAT7. *DNA Seq.* 16: 386-390.
7. Cao, J., Shan, D., Revett, T., Li, D., Wu, L., Liu, W., Tobin, J.F. and Gimeno, R.E. 2008. Molecular identification of a novel mammalian brain isoform of acyl-CoA:lysophospholipid acyltransferase with prominent ethanolamine lysophospholipid acylating activity, LPEAT2. *J. Biol. Chem.* 283: 19049-19057.
8. Online Mendelian Inheritance in Man, OMIM[™]. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 612039. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: Lpcat4 (mouse) mapping to 2 E3.

PRODUCT

LPAAT- η siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see LPAAT- η shRNA Plasmid (m): sc-149019-SH and LPAAT- η shRNA (m) Lentiviral Particles: sc-149019-V as alternate gene silencing products.

For independent verification of LPAAT- η (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-149019A, sc-149019B and sc-149019C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

LPAAT- η siRNA (m) is recommended for the inhibition of LPAAT- η expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor LPAAT- η gene expression knockdown using RT-PCR Primer: LPAAT- η (m)-PR: sc-149019-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.