

MS4A3 siRNA (m): sc-149647

BACKGROUND

MS4A (membrane-spanning 4-domain family, subfamily A) is a large family of proteins that includes at least 26 members in mouse and humans. Flanked by amino- and carboxyl-cytoplasmic regions, MS4A family members contain four highly conserved transmembrane domains. CD20, the most well-known MS4A family member, is a B-cell-specific molecule that functions as a calcium-permeable cation channel and is known to accelerate the G₀ to G₁ progression induced by IGF-1. MS4A3 (membrane-spanning four-domains, subfamily A, member 3), also known as HTM4 (hematopoietic-specific transmembrane protein 4) or CD20L, is a 214 amino acid protein that acts as a modulator for the G₁-S cell cycle transition. A multi-pass membrane protein, MS4A3 localizes to the endomembrane system and is also found in the cytoplasm. MS4A3 is expressed in hematopoietic cells and tissues and interacts with CDKN3.

REFERENCES

- Adra, C.N., Lelias, J.M., Kobayashi, H., Kaghad, M., Morrison, P., Rowley, J.D. and Lim, B. 1994. Cloning of the cDNA for a hematopoietic cell-specific protein related to CD20 and the β subunit of the high-affinity IgE receptor: evidence for a family of proteins with four membrane-spanning regions. *Proc. Natl. Acad. Sci. USA* 91: 10178-10182.
- Ishibashi, K., Suzuki, M., Sasaki, S. and Imai, M. 2001. Identification of a new multigene four-transmembrane family (MS4A) related to CD20, HTm4 and β subunit of the high-affinity IgE receptor. *Gene* 264: 87-93.
- Liang, Y. and Tedder, T.F. 2001. Identification of a CD20-, Fc ϵ RI β -, and HTm4-related gene family: sixteen new MS4A family members expressed in human and mouse. *Genomics* 72: 119-127.
- Liang, Y., Buckley, T.R., Tu, L., Langdon, S.D. and Tedder, T.F. 2001. Structural organization of the human MS4A gene cluster on chromosome 11q12. *Immunogenetics* 53: 357-368.
- Donato, J.L., Ko, J., Kutok, J.L., Cheng, T., Shirakawa, T., Mao, X.Q., Beach, D., Scadden, D.T., Sayegh, M.H. and Adra, C.N. 2002. Human HTm4 is a hematopoietic cell cycle regulator. *J. Clin. Invest.* 109: 51-58.
- Gray, C.H. and Barford, D. 2003. Getting in the ring: proline-directed substrate specificity in the cell cycle proteins Cdc14 and CDK2-cyclinA3. *Cell Cycle* 2: 500-502.
- Chinami, M., Yano, Y., Yang, X., Salahuddin, S., Moriyama, K., Shiroishi, M., Turner, H., Shirakawa, T. and Adra, C.N. 2005. Binding of HTm4 to cyclin-dependent kinase (Cdk)-associated phosphatase (KAP)-Cdk2-cyclin A complex enhances the phosphatase activity of KAP, dissociates cyclin A, and facilitates KAP dephosphorylation of Cdk2. *J. Biol. Chem.* 280: 17235-17242.

CHROMOSOMAL LOCATION

Genetic locus: Ms4a3 (mouse) mapping to 19 A.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

MS4A3 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see MS4A3 shRNA Plasmid (m): sc-149647-SH and MS4A3 shRNA (m) Lentiviral Particles: sc-149647-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

MS4A3 siRNA (m) is recommended for the inhibition of MS4A3 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor MS4A3 gene expression knockdown using RT-PCR Primer: MS4A3 (m)-PR: sc-149647-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.