



## NET-2 siRNA (m): sc-149913

### BACKGROUND

NET-2, also known as Tetraspanin-12 and Transmembrane 4 superfamily member 12, is a 305 amino acid multi-pass membrane protein that belongs to the transmembrane 4 superfamily, also known as the tetraspanin family. Members of the tetraspanin family are cell-surface proteins that are characterized by the presence of four hydrophobic domains and mediate signal transduction events that play a role in the regulation of cell development, activation, growth, motility, differentiation, and cancer. Considered molecular facilitators, tetraspanin proteins may regulate vesicle fusion and fission. Specifically, NET-2 plays a central role in retinal vascularization via regulation of NORRIN signal transduction. NET-2 also regulates membrane proteinases such as MT-MMP-1 and ADAM10. There are two isoforms of NET-2 which are produced as a result of alternative splicing events.

### REFERENCES

1. Todd, S.C., Doctor, V.S. and Levy, S. 1998. Sequences and expression of six new members of the tetraspanin/TM4SF family. *Biochim. Biophys. Acta* 1399: 101-104.
2. Serru, V., Dessen, P., Boucheix, C. and Rubinstein, E. 2000. Sequence and expression of seven new tetraspans. *Biochim. Biophys. Acta* 1478: 159-163.
3. Berditshevski, F. 2001. Complexes of tetraspanins with integrins: more than meets the eye. *J. Cell Sci.* 114: 4143-4151.
4. Hübner, K., Windoffer, R., Hutter, H. and Leube, R.E. 2002. Tetraspan vesicle membrane proteins: synthesis, subcellular localization, and functional properties. *Int. Rev. Cytol.* 214: 103-159.
5. Yunta, M. and Lazo, P.A. 2003. Tetraspanin proteins as organisers of membrane microdomains and signalling complexes. *Cell. Signal.* 15: 559-564.
6. Tarrant, J.M., Robb, L., van Spriel, A.B. and Wright, M.D. 2003. Tetraspanins: molecular organisers of the leukocyte surface. *Trends Immunol.* 24: 610-617.
7. Xu, D., Sharma, C. and Hemler, M.E. 2009. Tetraspanin-12 regulates ADAM10-dependent cleavage of amyloid precursor protein. *FASEB J.* 23: 3674-3681.
8. Lafleur, M.A., Xu, D. and Hemler, M.E. 2009. Tetraspanin proteins regulate membrane type-1 matrix metalloproteinase-dependent pericellular proteolysis. *Mol. Biol. Cell* 20: 2030-2040.
9. Online Mendelian Inheritance in Man, OMIM™. 2010. Johns Hopkins University, Baltimore, MD. MIM Number: 613138. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

### CHROMOSOMAL LOCATION

Genetic locus: Tspan12 (mouse) mapping to 6 A3.1.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

### PRODUCT

NET-2 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see NET-2 shRNA Plasmid (m): sc-149913-SH and NET-2 shRNA (m) Lentiviral Particles: sc-149913-V as alternate gene silencing products.

For independent verification of NET-2 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-149913A, sc-149913B and sc-149913C.

### STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

### APPLICATIONS

NET-2 siRNA (m) is recommended for the inhibition of NET-2 expression in mouse cells.

### SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

### RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor NET-2 gene expression knockdown using RT-PCR Primer: NET-2 (m)-PR: sc-149913-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.