

neugrin siRNA (m): sc-149923

BACKGROUND

Neugrin, also known as NGRN, mesenchymal stem cell protein DSC92, neurite outgrowth-associated protein or spinal cord-derived protein FI58G, is a 291 amino acid protein that plays a role in neuronal differentiation and belongs to the neugrin family. As both a secreted and nuclear protein, neugrin exists as two alternatively spliced isoforms and is highly expressed in skeletal muscle, brain and heart. Neugrin is upregulated in neuroblastoma cells by retinoic acid treatment and is encoded by a gene that maps to human chromosome 15q26.1. Chromosome 15 houses over 700 genes and comprises nearly 3% of the human genome. Angelman syndrome, Prader-Willi syndrome, Tay-Sachs disease and Marfan syndrome are all associated with defects in chromosome 15-localized genes.

REFERENCES

1. Roig Cayón, M., Basilio Monne, J. and Canalda Sahli, C. 1991. Manual instrumentation of root canals. Review of the last decade. Av. Odontostomatol. 7: 49-57.
2. Hurowitz, G.I., Silver, J.M., Brin, M.F., Williams, D.T. and Johnson, W.G. 1993. Neuropsychiatric aspects of adult-onset Tay-Sachs disease: two case reports with several new findings. J. Neuropsychiatry Clin. Neurosci. 5: 30-36.
3. Ishigaki, S., Niwa, J., Yoshihara, T., Mitsuma, N., Doyu, M. and Sobue, G. 2000. Two novel genes, human neugrin and mouse m-neugrin, are upregulated with neuronal differentiation in neuroblastoma cells. Biochem. Biophys. Res. Commun. 279: 526-533.
4. Midla, G.S. 2008. Diagnosis and management of patients with Marfan syndrome. JAAPA 21: 21-25.
5. Dan, B. 2009. Angelman syndrome: current understanding and research prospects. Epilepsia 50: 2331-2339.
6. Ferrer-Bolufer, I., Dalmau, J., Quiroga, R., Oltra, S., Orellana, C., Monfort, S., Roselló, M., De La Osa, A. and Martinez, F. 2009. Tyrosinemia type 1 and Angelman syndrome due to paternal uniparental isodisomy 15. J. Inherit. Metab. Dis. 32: S349-S353.
7. Wawrzik, M., Unmehopa, U.A., Swaab, D.F., van de Nes, J., Buiting, K. and Horsthemke, B. 2010. The C15orf2 gene in the Prader-Willi syndrome region is subject to genomic imprinting and positive selection. Neurogenetics 11: 153-161.

CHROMOSOMAL LOCATION

Genetic locus: Ngrn (mouse) mapping to 7 D3.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

neugrin siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see neugrin shRNA Plasmid (m): sc-149923-SH and neugrin shRNA (m) Lentiviral Particles: sc-149923-V as alternate gene silencing products.

For independent verification of neugrin (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-149923A, sc-149923B and sc-149923C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

neugrin siRNA (m) is recommended for the inhibition of neugrin expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor neugrin gene expression knockdown using RT-PCR Primer: neugrin (m)-PR: sc-149923-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.