

Nkx-6.3 siRNA (m): sc-150001

BACKGROUND

Members of the Nkx family of homeodomain proteins are key regulators of growth and development in several tissues, including brain, heart and pancreas. The Nkx-6 family is involved in the patterning of the pancreas and central nervous system and consists of three proteins: Nkx-6.1, Nkx-6.2 and Nkx-6.3. Nkx-6.1 is responsible for cellular differentiation in the ventral neural tube and spinal meninges in response to Shh. Nkx-6.2 is also expressed during neural tube development by neural progenitor cells. During development, Nkx-6.2 regulates interneuron fates by repressing the expression of Dbx1, a class I homeodomain transcription repressor. Nkx-6.3 is a 265 amino acid homeobox protein that shows selective expression in the duodenal and glandular endoderm, in contrast to Nkx-6.1 and Nkx-6.2 that are broadly expressed in the ventral positions of the developing CNS. Nkx-6.3 is required in differentiation of gastrin-producing G-cells in the stomach and antrum. There are two isoforms of Nkx-6.3 that are produced as a result of alternative splicing events.

REFERENCES

1. Pedersen, J.K., Nelson, S.B., Jorgensen, M.C., Henseleit, K.D., Fujitani, Y., Wright, C.V., Sander, M. and Serup, P. 2005. Endodermal expression of Nkx-6 genes depends differentially on Pdx1. *Dev. Biol.* 288: 487-501.
2. Nelson, S.B., Janiesch, C. and Sander, M. 2005. Expression of Nkx-6 genes in the hindbrain and gut of the developing mouse. *J. Histochem. Cytochem.* 53: 787-790.
3. Alanentalo, T., Chatonnet, F., Karlen, M., Sulniute, R., Ericson, J., Andersson, E. and Ahlgren, U. 2006. Cloning and analysis of Nkx-6.3 during CNS and gastrointestinal development. *Gene Expr. Patterns* 6: 162-170.
4. Zhao, S., Jiang, H., Wang, W. and Mao, B. 2007. Cloning and developmental expression of the *Xenopus* Nkx-6 genes. *Dev. Genes Evol.* 217: 477-483.
5. Hafner, B.P., Choi, M.Y., Shivdasani, R.A. and Rowitch, D.H. 2008. Expression and function of Nkx-6.3 in vertebrate hindbrain. *Brain Res.* 1222: 42-50.
6. Klink, R., Serup, P., Madsen, O.D. and Jørgensen, M.C. 2008. Specificity of four monoclonal anti-Nkx6-1 antibodies. *J. Histochem. Cytochem.* 56: 415-424.
7. Choi, M.Y., Romer, A.I., Wang, Y., Wu, M.P., Ito, S., Leiter, A.B. and Shivdasani, R.A. 2008. Requirement of the tissue-restricted homeodomain transcription factor Nkx-6.3 in differentiation of gastrin-producing G cells in the stomach antrum. *Mol. Cell. Biol.* 28: 3208-3218.
8. Online Mendelian Inheritance in Man, OMIM[™]. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 610772. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: Nkx6-3 (mouse) mapping to 8 A2.

PRODUCT

Nkx-6.3 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Nkx-6.3 shRNA Plasmid (m): sc-150001-SH and Nkx-6.3 shRNA (m) Lentiviral Particles: sc-150001-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Nkx-6.3 siRNA (m) is recommended for the inhibition of Nkx-6.3 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Nkx-6.3 gene expression knockdown using RT-PCR Primer: Nkx-6.3 (m)-PR: sc-150001-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.