

Papilin siRNA (m): sc-152016

BACKGROUND

Papilin, also known as PAPLN (Papilin, proteoglycan-like sulfated glycoprotein), is a 1,278 amino acid secreted protein that belongs to the Papilin family. Existing as six alternatively spliced isoforms, Papilin consists of a BPTI/Kunitz inhibitor domain, three Ig-like C2-type (immunoglobulin-like) domains, a PLAC domain and five TSP type-1 domains. Essential for embryonic development of *Drosophila melanogaster* and *Caenorhabditis elegans*, Papilin functions as an extracellular matrix glycoprotein that influences cell rearrangements and modulates metalloproteinases during organogenesis. Papilin interacts with several extracellular matrix components and ADAMTS enzymes. Papilin is encoded by a gene located on human chromosome 14q24.2. Chromosome 14 houses over 700 genes and comprises nearly 3.5% of the human genome.

REFERENCES

1. Campbell, A.G., Fessler, L.I., Salo, T. and Fessler, J.H. 1987. Papilin: a *Drosophila* proteoglycan-like sulfated glycoprotein from basement membranes. *J. Biol. Chem.* 262: 17605-17612.
2. Kramerova, I.A., Kawaguchi, N., Fessler, L.I., Nelson, R.E., Chen, Y., Kramerov, A.A., Kusche-Gullberg, M., Kramer, J.M., Ackley, B.D., Sieron, A.L., Prockop, D.J. and Fessler, J.H. 2000. Papilin in development; a pericellular protein with a homology to the ADAMTS metalloproteinases. *Development* 127: 5475-5485.
3. Kramerova, I.A., Kramerov, A.A. and Fessler, J.H. 2003. Alternative splicing of Papilin and the diversity of *Drosophila* extracellular matrix during embryonic morphogenesis. *Dev. Dyn.* 226: 634-642.
4. Tucker, R.P. 2004. The thrombospondin type 1 repeat superfamily. *Int. J. Biochem. Cell Biol.* 36: 969-974.
5. Fessler, J.H., Kramerova, I., Kramerov, A., Chen, Y. and Fessler, L.I. 2004. Papilin, a novel component of basement membranes, in relation to ADAMTS metalloproteinases and ECM development. *Int. J. Biochem. Cell Biol.* 36: 1079-1084.
6. Kimura, K., Wakamatsu, A., Suzuki, Y., Ota, T., Nishikawa, T., Yamashita, R., Yamamoto, J., Sekine, M., Tsuritani, K., Wakaguri, H., Ishii, S., Sugiyama, T., Saito, K., Isono, Y., Irie, R., Kushida, N., Yoneyama, T., Otsuka, R., et al. 2006. Diversification of transcriptional modulation: large-scale identification and characterization of putative alternative promoters of human genes. *Genome Res.* 16: 55-65.
7. Kawano, T., Zheng, H., Merz, D.C., Kohara, Y., Tamai, K.K., Nishiwaki, K. and Culotti, J.G. 2009. *C. elegans* mig-6 encodes Papilin isoforms that affect distinct aspects of DTC migration, and interacts genetically with mig-17 and collagen IV. *Development* 136: 1433-1442.

CHROMOSOMAL LOCATION

Genetic locus: Papln (mouse) mapping to 12 D1.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Papilin siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Papilin shRNA Plasmid (m): sc-152016-SH and Papilin shRNA (m) Lentiviral Particles: sc-152016-V as alternate gene silencing products.

For independent verification of Papilin (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152016A, sc-152016B and sc-152016C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Papilin siRNA (m) is recommended for the inhibition of Papilin expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Papilin gene expression knockdown using RT-PCR Primer: Papilin (m)-PR: sc-152016-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.