

PCPE-2 siRNA (m): sc-152117

BACKGROUND

PCPE-2 (procollagen C-endopeptidase enhancer 2), also known as procollagen COOH-terminal proteinase enhancer 2, is a 415 amino acid secreted protein that binds to the C-terminal propeptides of type I and II procollagens. PCPE-2 contains two CUB domains, which are involved in protein-protein interactions and an NTR domain, which potentially is involved in the regulation of metalloproteases such as BMP-1. Due to its strong affinity to heparin, PCPE-2 may play a role in the blood coagulation system. PCPE-2 is expressed in heart, lungs, brain, kidney, skeletal muscle, placenta, pancreas and trabecular meshwork. Interestingly, PCPE-2 enhances cleavage of the hexapeptide extension of pro-apolipoprotein (apo)-AI, thereby accelerating the proteolytic processing of apoAI. This evidence suggests that PCPE-2 may be involved in the regulation of apoAI synthesis and HDL levels, potentially playing a role in cardiovascular disease.

REFERENCES

1. Srivastava, R.A. 1994. Saturated fatty acid, but not cholesterol, regulates apolipoprotein AI gene expression by posttranscriptional mechanism. *Biochem. Mol. Biol. Int.* 34: 393-402.
2. Takahara, K., Osborne, L., Elliott, R.W., Tsui, L.C., Scherer, S.W. and Greenspan, D.S. 1996. Fine mapping of the human and mouse genes for the type I procollagen COOH-terminal proteinase enhancer protein. *Genomics* 31: 253-256.
3. Scott, I.C., Clark, T.G., Takahara, K., Hoffman, G.G. and Greenspan, D.S. 1999. Structural organization and expression patterns of the human and mouse genes for the type I procollagen COOH-terminal proteinase enhancer protein. *Genomics* 55: 229-234.
4. Xu, H., Acott, T.S. and Wirtz, M.K. 2000. Identification and expression of a novel type I procollagen C-proteinase enhancer protein gene from the glaucoma candidate region on 3q21-q24. *Genomics* 66: 264-273.
5. Steiglit, B.M., Keene, D.R. and Greenspan, D.S. 2002. PCOLCE2 encodes a functional procollagen C-proteinase enhancer (PCPE2) that is a collagen-binding protein differing in distribution of expression and post-translational modification from the previously described PCPE1. *J. Biol. Chem.* 277: 49820-49830.
6. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607064. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. Hartigan, N., Garrigue-Antar, L. and Kadler, K.E. 2003. Bone morphogenetic protein-1 (BMP-1). Identification of the minimal domain structure for procollagen C-proteinase activity. *J. Biol. Chem.* 278: 18045-18049.
8. Zhu, J., Gardner, J., Pullinger, C.R., Kane, J.P., Thompson, J.F. and Franccone, O.L. 2009. Regulation of apoAI processing by procollagen C-proteinase enhancer-2 and bone morphogenetic protein-1. *J. Lipid Res.* 50: 1330-1339.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: Pcolce2 (mouse) mapping to 9 E3.3.

PRODUCT

PCPE-2 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see PCPE-2 shRNA Plasmid (m): sc-152117-SH and PCPE-2 shRNA (m) Lentiviral Particles: sc-152117-V as alternate gene silencing products.

For independent verification of PCPE-2 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152117A, sc-152117B and sc-152117C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

PCPE-2 siRNA (m) is recommended for the inhibition of PCPE-2 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor PCPE-2 gene expression knockdown using RT-PCR Primer: PCPE-2 (m)-PR: sc-152117-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.