

pierce 1 siRNA (m): sc-152249

BACKGROUND

Pierce 1, also known as UPF0691 protein C9orf116, is a 136 amino acid protein that belongs to the UPF0691 family. Pierce 1 exists as two alternatively spliced isoforms, which are conserved in chimpanzee, canine, bovine, mouse and rat. Pierce 1 is predominantly expressed in micrometastases and macrometastases of human small cell lung cancer. Pierce 1 is encoded by a gene that maps to chromosome 9q34.3, which consists of about 145 million bases, makes up 4% of the human genome and encodes nearly 900 genes. Notably, chromosome 9 encompasses the largest interferon family gene cluster. Considered to play a role in gender determination, deletion of the distal portion of 9p can lead to development of male to female sex reversal, the phenotype of a female with a male X,Y genotype. Hereditary hemorrhagic telangiectasia, which is characterized by harmful vascular defects, is associated with chromosome 9.

REFERENCES

1. Wjst, M., Fischer, G., Immervoll, T., Jung, M., Saar, K., Rueschendorf, F., Reis, A., Ulbrecht, M., Gomolka, M., Weiss, E.H., Jaeger, L., Nickel, R., Richter, K., Kjellman, N.I., Griese, M., von Berg, A., Gappa, M., et al. 1999. A genome-wide search for linkage to asthma. German Asthma Genetics Group. *Genomics* 58: 1-8.
2. Kakiuchi, S., Daigo, Y., Tsunoda, T., Yano, S., Sone, S. and Nakamura, Y. 2003. Genome-wide analysis of organ-preferential metastasis of human small cell lung cancer in mice. *Mol. Cancer Res.* 1: 485-499.
3. Humphray, S.J., Oliver, K., Hunt, A.R., Plumb, R.W., Loveland, J.E., Howe, K.L., Andrews, T.D., Searle, S., Hunt, S.E., Scott, C.E., Jones, M.C., Ainscough, R., Almeida, J.P., Ambrose, K.D., Ashwell, R.I., et al. 2004. DNA sequence and analysis of human chromosome 9. *Nature* 429: 369-374.
4. Fernandez-L, A., Garrido-Martin, E.M., Sanz-Rodriguez, F., Pericacho, M., Rodriguez-Barbero, A., Eleno, N., Lopez-Novoa, J.M., Düwell, A., Vega, M.A., Bernabeu, C. and Botella, L.M. 2007. Gene expression fingerprinting for human hereditary hemorrhagic telangiectasia. *Hum. Mol. Genet.* 16: 1515-1533.
5. Cottin, V., Dupuis-Girod, S., Lesca, G. and Cordier, J.F. 2007. Pulmonary vascular manifestations of hereditary hemorrhagic telangiectasia (Rendu-Osler disease). *Respiration* 74: 361-378.
6. Baines, K.J., Wood, L.G. and Gibson, P.G. 2009. The nutrigenomics of asthma: molecular mechanisms of airway neutrophilia following dietary antioxidant withdrawal. *OMICS* 13: 355-365.
7. SWISS-PROT/TrEMBL (Q5BN46). World Wide Web URL: <http://www.uniprot.org/uniprot/Q5BN46>

CHROMOSOMAL LOCATION

Genetic locus: 1700007K13Rik (mouse) mapping to 2 A3.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

pierce 1 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see pierce 1 shRNA Plasmid (m): sc-152249-SH and pierce 1 shRNA (m) Lentiviral Particles: sc-152249-V as alternate gene silencing products.

For independent verification of pierce 1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152249A and sc-152249B.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

pierce 1 siRNA (m) is recommended for the inhibition of pierce 1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor pierce 1 gene expression knockdown using RT-PCR Primer: pierce 1 (m)-PR: sc-152249-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.