



Polycystin-1L2 siRNA (m): sc-152384

BACKGROUND

Polycystin-1L2, also known as PKD1L2 or PC1L2, is a 2,459 amino acid multi-pass membrane protein that belongs to the Polycystin family. Expressed in a variety of tissues, including placenta, brain, liver, lung, testis, skeletal muscle and fetal and adult heart, Polycystin-1L2 is thought to function as an ion-channel regulator and may also exhibit activity as a G protein-coupled receptor. Polycystin-1L2 contains a latrophilin/CL-1-like GPCR proteolytic site (GPS) domain, a Polycystin-1, lipoxigenase, α -toxin (PLAT) domain and several transmembrane domains through which it conveys its regulatory function. Human Polycystin-1L2 shares 73% sequence similarity with its mouse counterpart, suggesting a conserved role between species. Defects in the gene encoding Polycystin-1L2 may be associated with Polycystic kidney disease, a progressive disorder characterized by the presence of cysts in the kidneys. Seven isoforms of Polycystin-1L2 exist due to alternative splicing events.

REFERENCES

1. Wu, G., Hayashi, T., Park, J.H., Dixit, M., Reynolds, D.M., Li, L., Maeda, Y., Cai, Y., Coca-Prados, M. and Somlo, S. 1998. Identification of PKD2L, a human PKD2-related gene: tissue-specific expression and mapping to chromosome 10q25. *Genomics* 54: 564-568.
2. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607894. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Li, A., Tian, X., Sung, S.W. and Somlo, S. 2003. Identification of two novel polycystic kidney disease-1-like genes in human and mouse genomes. *Genomics* 81: 596-608.
4. Yuasa, T., Takakura, A., Denker, B.M., Venugopal, B. and Zhou, J. 2004. Polycystin-1L2 is a novel G-protein-binding protein. *Genomics* 84: 126-138.
5. Augustyniak-Bartosik, H., Weyde, W., Krajewska, M. and Klinger, M. 2004. Autosomal dominant polycystic kidney disease—research status and clinical manifestation. *Postepy Hig. Med. Dosw.* 58: 530-537.
6. Hovater, M.B., Olteanu, D., Welty, E.A. and Schwiebert, E.M. 2008. Purinergic signaling in the lumen of a normal nephron and in remodeled PKD encapsulated cysts. *Purinergic Signal*. 4: 109-124.
7. Aguiari, G., Trimi, V., Bogo, M., Mangolini, A., Szabadkai, G., Pinton, P., Witzgall, R., Harris, P.C., Borea, P.A., Rizzuto, R. and Del Senno, L. 2008. Novel role for Polycystin-1 in modulating cell proliferation through calcium oscillations in kidney cells. *Cell Prolif.* 41: 554-573.

CHROMOSOMAL LOCATION

Genetic locus: Pkd1l2 (mouse) mapping to 8 E1.

PRODUCT

Polycystin-1L2 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Polycystin-1L2 shRNA Plasmid (m): sc-152384-SH and Polycystin-1L2 shRNA (m) Lentiviral Particles: sc-152384-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Polycystin-1L2 siRNA (m) is recommended for the inhibition of Polycystin-1L2 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Polycystin-1L2 gene expression knockdown using RT-PCR Primer: Polycystin-1L2 (m)-PR: sc-152384-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.