



Polycystin-1L3 siRNA (m): sc-152385

BACKGROUND

Polycystin-1L3, also known as PKD1L3 (polycystic kidney disease 1-like 3), is a 1,732 amino acid multi-pass membrane protein that contains one PLAT domain, one GPS domain and one C-type lectin domain. Expressed at high levels in placenta and present at lower levels in lung and heart, Polycystin-1L3 is thought to function as an ion-channel regulator that may interact with Polycystin-L and play a role in heteromeric taste channels. The gene encoding Polycystin-1L3 maps to human chromosome 16, which encodes over 900 genes and comprises nearly 3% of the human genome. The GAN gene is located on chromosome 16 and, with mutation, may lead to giant axonal neuropathy, a nervous system disorder characterized by increasing malfunction with growth. The rare disorder Rubinstein-Taybi syndrome is also associated with chromosome 16, as is Crohn's disease, which is a gastrointestinal inflammatory condition.

REFERENCES

- Li, A., Tian, X., Sung, S.W. and Somlo, S. 2003. Identification of two novel polycystic kidney disease-1-like genes in human and mouse genomes. *Genomics* 81: 596-608.
- LopezJimenez, N.D., Cavenagh, M.M., Sainz, E., Cruz-Ithier, M.A., Battey, J.F. and Sullivan, S.L. 2006. Two members of the TRPP family of ion channels, Pkd1l3 and Pkd2l1, are co-expressed in a subset of taste receptor cells. *J. Neurochem.* 98: 68-77.
- Ishimaru, Y., Inada, H., Kubota, M., Zhuang, H., Tominaga, M. and Matsu-nami, H. 2006. Transient receptor potential family members PKD1L3 and PKD2L1 form a candidate sour taste receptor. *Proc. Natl. Acad. Sci. USA* 103: 12569-12574.
- Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 607895. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Chen, W.C., Tzeng, Y.S. and Li, H. 2008. Gene expression in early and progression phases of autosomal dominant polycystic kidney disease. *BMC Res Notes* 1: 131.

CHROMOSOMAL LOCATION

Genetic locus: Pkd1l3 (mouse) mapping to 8 D3.

PRODUCT

Polycystin-1L3 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Polycystin-1L3 shRNA Plasmid (m): sc-152385-SH and Polycystin-1L3 shRNA (m) Lentiviral Particles: sc-152385-V as alternate gene silencing products.

For independent verification of Polycystin-1L3 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152385A, sc-152385B and sc-152385C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Polycystin-1L3 siRNA (m) is recommended for the inhibition of Polycystin-1L3 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Polycystin-1L3 gene expression knockdown using RT-PCR Primer: Polycystin-1L3 (m)-PR: sc-152385-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.