

Rab11-FIP5 siRNA (m): sc-152664

BACKGROUND

RAB11 is a GTPase that regulates endosomal trafficking to apical plasma membrane domains in polarized epithelial cells. Rab11-FIP5 (RAB11 family interacting protein 5 (class I)) is a 653 amino acid cytoplasmic protein that is a Rab effector involved in protein trafficking from apical recycling endosomes to the apical plasma membrane. Rab11-FIP5 forms an heterooligomeric complex with Rab11-FIP4, binds γ -SNAP and 60 kDa Ro/SSA, and binds Rab 11A that has been activated by GTP binding. Rab11-FIP5 is detected at low levels in heart, brain, placenta, lung, liver, adipocytes, kidney, spleen, skeletal muscle and pancreas, and contains one C2 domain and one FIP-RBD domain. The Rab11-FIP5 gene is conserved in canine, bovine, mouse, chicken and zebrafish, and maps to human chromosome 2p13.2.

REFERENCES

1. Wang, D., Buyon, J.P., Zhu, W. and Chan, E.K. 1999. Defining a novel 75-kDa phosphoprotein associated with SS-A/Ro and identification of distinct human autoantibodies. *J. Clin. Invest.* 104: 1265-1275.
2. Prekeris, R., Klumperman, J. and Scheller, R.H. 2000. A Rab11/Rip11 protein complex regulates apical membrane trafficking via recycling endosomes. *Mol. Cell* 6: 1437-1448.
3. Chen, D., Xu, W., He, P., Medrano, E.E. and Whiteheart, S.W. 2001. Gaf-1, a γ -SNAP-binding protein associated with the mitochondria. *J. Biol. Chem.* 276: 13127-13135.
4. Prekeris, R., Davies, J.M. and Scheller, R.H. 2001. Identification of a novel Rab11/25 binding domain present in Eferin and Rip proteins. *J. Biol. Chem.* 276: 38966-38970.
5. Hales, C.M., Griner, R., Hobdy-Henderson, K.C., Dorn, M.C., Hardy, D., Kumar, R., Navarre, J., Chan, E.K., Lapierre, L.A. and Goldenring, J.R. 2001. Identification and characterization of a family of Rab11-interacting proteins. *J. Biol. Chem.* 276: 39067-39075.
6. Online Mendelian Inheritance in Man, OMIM[™]. 2001. Johns Hopkins University, Baltimore, MD. MIM Number: 605536. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. Roohi, J., Tegay, D.H., Pomeroy, J.C., Burkett, S., Stone, G., Stanyon, R. and Hatchwell, E. 2008. A *de novo* apparently balanced translocation [46,XY,t(2;9)(p13;p24)] interrupting RAB11FIP5 identifies a potential candidate gene for autism spectrum disorder. *Am. J. Med. Genet. B Neuropsychiatr. Genet.* 147B: 411-417.
8. Schonteich, E., Wilson, G.M., Burden, J., Hopkins, C.R., Anderson, K., Goldenring, J.R. and Prekeris, R. 2008. The Rip11/Rab11-FIP5 and kinesin II complex regulates endocytic protein recycling. *J. Cell Sci.* 121: 3824-3833.

CHROMOSOMAL LOCATION

Genetic locus: Rab11fip5 (mouse) mapping to 6 C3.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Rab11-FIP5 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Rab11-FIP5 shRNA Plasmid (m): sc-152664-SH and Rab11-FIP5 shRNA (m) Lentiviral Particles: sc-152664-V as alternate gene silencing products.

For independent verification of Rab11-FIP5 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152664A, sc-152664B and sc-152664C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Rab11-FIP5 siRNA (m) is recommended for the inhibition of Rab11-FIP5 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Rab11-FIP5 gene expression knockdown using RT-PCR Primer: Rab11-FIP5 (m)-PR: sc-152664-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.